
# Invacare® Soft Tilt

en	Accessory for Invacare® SB® 755 User Manual .....	3
da	Tilbehør til Invacare® SB® 755 Brugsanvisning .....	25
no	Tilbehør for Invacare® SB® 755 Bruksanvisning .....	49
sv	Tillbehör för Invacare® SB® 755 Bruksanvisning .....	71

This manual MUST be given to the user of the product.  
BEFORE using this product, read this manual and save for future  
reference.

**Yes, you can.®**

©2017 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

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# Contents

---

This manual MUST be given to the user of the product.  
BEFORE using this product, read this manual and save for future reference.

<b>1 General</b> . . . . .	<b>4</b>
1.1 Introduction . . . . .	4
1.1.1 Definitions . . . . .	4
1.2 Symbols in this manual . . . . .	4
1.3 Intended use. . . . .	4
1.4 Service life . . . . .	4
1.5 Warranty information . . . . .	4
1.6 Limitation of liability . . . . .	4
<b>2 Safety</b> . . . . .	<b>5</b>
2.1 General safety information . . . . .	5
2.2 Labels and symbols on the product . . . . .	6
2.2.1 Product label . . . . .	6
<b>3 Setup</b> . . . . .	<b>7</b>
3.1 General Safety Information . . . . .	7
3.2 Scope of delivery . . . . .	7
3.3 Mounting the Soft Tilt . . . . .	7
<b>4 Usage</b> . . . . .	<b>10</b>
4.1 General Safety Information . . . . .	10
4.2 Overview . . . . .	10
4.3 Control Box and Hand Control . . . . .	10
4.4 On/Off and CPR/Reset function . . . . .	10
4.5 Emergency Stop . . . . .	11
4.6 Manual operation . . . . .	11
4.6.1 Recommended Procedures . . . . .	11
4.7 Automatic operation . . . . .	11
4.7.1 Auto program sequence . . . . .	12
<b>5 Programming</b> . . . . .	<b>13</b>
5.1 Individually defined programs . . . . .	13
5.2 How to program the individual automated programs . . . . .	14
5.3 Program examples for automated pressure relief . . . . .	15
<b>6 Maintenance</b> . . . . .	<b>17</b>
6.1 General Maintenance Information . . . . .	17
6.2 Cleaning and disinfection. . . . .	17
6.3 Battery . . . . .	17
<b>7 After Use</b> . . . . .	<b>18</b>
7.1 Waste disposal . . . . .	18
<b>8 Troubleshooting</b> . . . . .	<b>19</b>
8.1 Troubleshooting the electrical system. . . . .	19
8.2 Programming Errors. . . . .	20
<b>9 Technical Data</b> . . . . .	<b>21</b>
9.1 Characteristics. . . . .	21
9.2 Weights . . . . .	21
9.3 Allowed mattress sizes . . . . .	21
9.4 Electrical system . . . . .	21
<b>10 Electromagnetic compatibility (EMC)</b> . . . . .	<b>22</b>
10.1 General EMC information . . . . .	22
10.2 Electromagnetic emission. . . . .	22
10.3 Electromagnetic Immunity . . . . .	22

# 1 General

## 1.1 Introduction

This manual belongs to an Invacare accessory and it contains important information about handling and assembly. To ensure safety when using the product, read the user manual of the related product carefully and follow the safety instructions.

Find the user manual on Invacare's website or contact your Invacare representative (see addresses at the end of this manual).

### 1.1.1 Definitions

All references to left and right are based on a person lying on his back in the bed, with his head in the head end.

## 1.2 Symbols in this manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

#### WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

#### CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

#### IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

#### Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

## 1.3 Intended use

The Soft Tilt and control unit is intended to be used only in conjunction with 85 or 90 cm wide Invacare® SB 755 beds, equipped with side rails in an appropriate height (see chart in *9 Technical Data, page 21*), side rail covers and castellated foam mattress (see size chart *9 Technical Data, page 21*)

- As a part of an overall pressure ulcer prevention programme of care.
- To move bedridden end-users to lay on left or right side or on the back.
- Support the end-user in getting in and out of bed.
- When used in the "automatic" mode it has been designed to provide pressure reduction to end-users.
- When used in the "manual" mode it is designed to support the care staff in turning the end-user from side to side in the bed.
- It is suitable for use in all home care, residential and nursing care settings.

## Indications

The intended patient group for Soft Tilt is:

- Bedridden immobile/partly immobile adult end-users.
- Who needs pressure care.
- Who needs to be nursed in the care bed such as personal hygiene, dressing etc. This group of end users will often be transferred from bed to chair using a lifter or manually supported by care staff in combination with a technical aid such as a rollator.

## Contraindications

The Soft Tilt is not intended for psychiatric, restless or anxious end-users or persons with spasm.

## Precautions

- Before Soft Tilt is used it is important to evaluate, if Soft Tilt is suitable for the actual end-user.
- For security reasons the end-user must not operate Soft Tilt themselves.

#### WARNING!

Any other or incorrect use could lead to hazardous situations. Invacare accepts no liability for any use, change or assembly of the product, other than stated in this user manual.

## 1.4 Service life

The expected service life of this product is five years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

## 1.5 Warranty information

Terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

Contact information for your local Invacare office is located at the end of this manual.

## 1.6 Limitation of liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorized modifications and/or use of unsuitable spare parts

## 2 Safety

### 2.1 General safety information

#### WARNING!

##### Compatibility with Invacare® SB 755 beds only

The Soft Tilt is specially designed and must only be used in conjunction with 85 or 90 cm wide Invacare® SB 755 beds.

#### WARNING!

##### Risk of falling

If the Soft Tilt is used without or with the wrong side rails, there is a risk for the end-user to fall out of the bed.

- Never use the Soft Tilt on a bed without side rails.
- Only use the Soft Tilt in combination with Invacare side rails listed in *9 Technical Data, page 21*.
- Always make sure the side rail is up in the opposite site of where the care stuff is placed OR one care staff is placed on each side of the bed, when the manual mode is used.
- Always make sure both side rails are up, when the automatic mode is activated.

#### WARNING!

##### Risk of squeezing and/or suffocation

If the side rails are used without a cover or with a non-breathable cover, there is a risk of squeezing and/or suffocation for the end-user.

- Always use a breathable cover for the side rails when using the Soft Tilt.

#### WARNING!

##### Risk of falling and/or squeezing

- Always perform a risk assessment of the end-user's condition and ability of moving.
- Never use the Soft Tilt if the end-user is anxious or restless.
- Always make sure the end-user is placed in the middle of the bed and all body parts are on the mattress before the wings are moved.
- Never leave the end-user unattended when the manual mode is used.
- The end-user must never get in or out of the bed while the wings are moved manually or the automatic mode is activated.
- The hand control must always be used by care giver.

#### WARNING!

##### Risk of injury or damage to property

Inappropriate handling of cables can cause electrical shock and product failure.

- Do not kink, shear or otherwise damage the main power cord.
- Do not roll the castors over the main power cord.
- Do not bring main power cord into moving parts.
- Disconnect the plug from the mains before moving the bed.
- Make sure that no cables (mains or from other equipment) are jammed or damaged, when the bed is used.
- Keep bed components and accessories at least 30 cm away from a heated surface and not in direct sunlight.

### Mattresses

#### WARNING!

##### Safety aspects regarding combination of side rails and mattresses:

To get the highest possible safety level, when using side rails on the bed, the minimum and maximum measures for mattresses, must be respected.


- For correct mattress measures see mattress tables in chapter *9 Technical Data, page 21*.

#### WARNING!

##### Risk of entrapment and/or suffocation

- The end-user could get trapped and/or suffocate, if the horizontal space, between the mattress side and the inside of the side rail, is too big. Follow the minimum width (and length) of mattresses in combination with a side rail, as stated in the mattress table in chapter *9 Technical Data, page 21*.
- Be aware that using very thick or soft mattresses (low density), or a combination of these, increases the risk.

#### WARNING!

##### Risk of falling

The end-user can fall over the edge and get seriously injured, if the vertical distance  $\textcircled{A}$  between the top of the mattress and the edge of the side rail/bed end, is too short. See image above.

- Always keep a minimum distance  $\textcircled{A}$  of 22 cm on the side of the wing the user is lying on (secondary wing elevated to max. 12°).
- Follow the maximum mattress height in combination with the side rail as stated in the mattress table in chapter *9 Technical Data, page 21*.

## Electromagnetic interference

### WARNING!

#### Risk of malfunction due to electromagnetic interference

Electromagnetic interference between this product and other electrical equipment can occur and disturb the electrical adjustment functions of this product. To prevent, reduce or eliminate such electromagnetic interference:

- Only use original cables, accessories and spare parts, to not increase electromagnetic emission or reduce electromagnetic immunity of this product.
- Do not use portable radio frequency (RF) communications equipment closer than 30 cm to any part of this product (including cables).
- Do not use this product near active high-frequency surgical equipment and the RF shielded room of a system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.
- If disturbances occur, increase the distance between this product and the other equipment or switch it off.
- Refer to the detailed information and follow the guidance in chapter 10 *Electromagnetic compatibility (EMC)*, page 22.


The product label is placed on the main module of the Soft Tilt and contains the main product information, including technical data.

### Symbols

	Serial Number
	Reference Number
	Manufacturer Address
	Manufacturing Date
	Max. User Weight
	Max. Safe Working Load
	Insulation Class 2
	Type B Applied Part
	WEEE conform
	This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

## 2.2 Labels and symbols on the product

### 2.2.1 Product label

### Soft Tilt for Invacare® SB 755

#### Abbreviations for technical data:

- Iin = Incoming Current
- AC = Alternating Current
- Uin = Incoming Voltage
- Max = maximum
- Int. = Intermittence
- min = minutes

For more information about technical data, refer to 9 *Technical Data*, page 21.

## 3 Setup

### 3.1 General Safety Information

When you receive the product, check the packaging. If the packaging shows any signs of damage upon delivery, contact the shipping company.

#### WARNING!

#### Risk of injury or damage to property


The assembly and installation must be done by authorized or trained personnel.

- Follow instructions carefully. If you have any question regarding the assembly, contact your local dealer or Invacare representative.
- If the product shows any sign of damages, do not use the product. Contact your local dealer or Invacare representative.
- The electrical equipment of the product must not be dismantled or combined with other electrical equipment.
- After each assembly, check that all fittings are properly tightened and that all parts have the correct function.

### 3.2 Scope of delivery

The Soft Tilt is delivered partially assembled in a cardboard box.

#### Main Parts


- Ⓐ Head end arm of left wing
  - Ⓑ Head end arm of right wing
  - Ⓒ Main module (pre-assembled unit)
  - Ⓓ Foot end arm of left wing
  - Ⓔ Foot end arm of right wing
  - Ⓕ 2 sets: Fixation plate, 2x carriage-bolt, 2x washer, 2x wing-nut
  - Ⓖ 4 sets: Allen-bolt, washer, wing-nut
  - Ⓗ Control Box
  - Ⓘ Hand Control
- Textile cover (not show in image)

### 3.3 Mounting the Soft Tilt

#### Mounting the main module

1.


a.

On each side of the bed, remove the nut Ⓐ and bolt Ⓑ using two 13 mm keys to remove the slats assembly Ⓒ from the bed.


- ⓘ For beds equipped with wooden slats, 2 spare aluminium slats are required to replace the wooden slats of the slats assembly Ⓒ.

b. Additional step only required for dividable beds (SB 755D)


Unscrew and remove the distance stick, located at the cross bar of the mattress support head section, and keep it available for remounting, when the Soft Tilt is removed again from the bed.

2.

Place the four carriage bolts Ⓜ into the holes of the brackets on the main module Ⓒ and remove one of the brackets Ⓝ from the slats assembly Ⓒ.


3.

Insert the loose ends of the two slats underneath the cover plates of the main module and insert them back into the bracket N on the other side.


4.

**! IMPORTANT!**  
Do not hold the main module on the middle plate while lifting, as it can become loose.

Place the main module with slats assembly onto the crossbars P, with the long connection arms © facing towards the head end of the bed. Ensure that the cables of the motors run between the crossbars P and no cables are squeezed. Fix the slats assembly with bolts K and nuts J removed in step 1.

5.

Ensure that the main module is placed in the middle of the mattress support frame and the distance to the side tubes is the same on both sides. Fix the main module to the crossbars of the mattress support frame, with the fixation plate washers and nuts, from underneath on both sides of the bed.

### Wiring

The control box and hand control is to be hung at the footboard.

#### Connections at the control box

- © Power supply
- ® Hand control and squeeze protection switches via splitter cable
- © Right motor
- © Left motor

1. Connect the cables of the right and left motor to the two extension cables © and © from the control box.
2. Connect the hand control to the shorter end of the splitter cable ® and connect the power cable © to the mains.
3. Deactivate the emergency stop, turn the system on, activate manual mode and use the hand control to check that the left and right motor are correctly connected to the control box (see chapter *Usage* for detailed instructions).


4.

Elevate both wings of the main module and remove the tape fixating the cable of the squeeze protection switches on the main module. Run the switch cables to one side of the bed and run the free end underneath the first crosstube and on top of the second crosstube towards the foot end of the bed. Connect the cable to the longest end of the splitter cable ®. Fix each of the cables from the squeeze protection switches with a plastic strip at the crosstube and Soft Tilt frame, facing the head end of the bed.

5.


Ensure that all cables run underneath the bed slats and fasten one motor cable together with the switch cable on side of the bed and the second motor cable on the other side of the bed, with three plastic strips on each side. The first strip is placed on the crossbar and the remaining two on the side tube towards the foot end of the bed.

 Move the bed all the way up and down to ensure no cables are getting squeezed.

### Mounting the arms and cover


1. Ensure both wings of the main module are slightly elevated, so the connections for the arms are accessible.
- 2.

Connect the four arms with the main module, ensuring the arms with the hinge is placed over the knee bend hinge on the bed and that the curved shaped sides of all arms are pointing outwards. Fix the four arms with the Allen bolts and tighten the wing-nuts facing inwards. Ensure the safety washers are placed underneath the wing-nuts.


3.

**!** **IMPORTANT!**  
When mounting, ensure that the four arms are in a horizontal position and adjust them if necessary

To adjust the position of the arms, fold them up at the joints to the main module and adjust screw  until the arms are in a horizontal position when folded back down.

4.

Mount the cover by sliding the pockets over the head end and foot end arms. The head end of the cover is marked by "TOP".

**!** **IMPORTANT!**  
The Soft Tilt must not be used without the cover.

5. Place and secure the mattress onto the Soft Tilt.
  - a. Release the Velcro fastener of both straps on the side of the cover and fold them away to the side.
  - b. Place the mattress onto the bed on top of the Soft Tilt (for compatible mattresses see 9 *Technical Data*, page 21 ).
  - c. Guide the straps over the mattress and re-fix the Velcro fastener to tightly secure the mattress to the Soft Tilt.

## 4 Usage

### 4.1 General Safety Information

**WARNING!**

**Risk of personal injury and damage to property.**

- The bed must be placed so that the height adjustment is not obstructed by, for example, lifts or furniture.
- Take care that no body parts are being squeezed between fixed parts (such as side rails, bed ends etc) and moving parts.
- The hand control must not be used by children.
- The hand control must only be used by care giver.

**IMPORTANT**

If a power failure occurs, the battery secures the possibility of resetting the Soft Tilt via the “CPR” function.

- Make sure the battery is fully charged before using the Soft Tilt.
- refer to Battery in the Maintenance section for further information.

### 4.2 Overview

The Soft Tilt can be operated either in manual or automatic mode.

In both modes it is possible to activate the normal functions the bed.

It is recommend to have the head section slightly raised, while using the Soft Tilt, to provide a more comfortable position to the end-user.

#### Manual mode

The manual mode has been designed to support the care staff in turning the end-user from side to side in the bed.

The wings can be raised and lowered via the hand control (see section *Manual Operation*).


#### Automatic mode

The automatic mode has been designed to provide pressure reduction to end-users.

In the AUTO program, the wings will move automatically according to pre-defined parameters for time and angle (see section *Auto Program Sequence*).

In Program 1 and 2, the wings will move automatically according to individually defined parameters for time and angle (see chapter *Programming*).

### 4.3 Control Box and Hand Control

- (A) On/Off & CPR (cardiac arrest function) button
- (B) Emergency stop button
- (C) Function buttons to choose the mode of operation.
- (D) Setup buttons to set the parameters for program 1 and 2
- (E) Display
- (F) Position indicators
- (G) Hand control

### 4.4 On/Off and CPR/Reset function

#### Turn the system On

1. Press and hold the button (A) for 3 sec. to turn the system on.

#### CPR / Reset function

The CPR (cardiopulmonary resuscitation) / Reset function overrides all commands and moves both wings into a horizontal position.

1. To activate the CPR / Reset function press button (A) (without holding).

#### Turn the system off

1. Press the button (A) to reset the system and if applicable, wait until both wings moved into horizontal position.
2. Press and hold the button (A) for 3 sec. to turn the system off.

#### Display

System Status	Display	Note
OFF		Display blank
ON	software version (3 digit number)	Displayed 2 sec. after system has been turned on. (switches to passive mode if no function is selected).
ON	— — —	Passive mode (no function selected)
Reset	OFF	Displayed after the CPR / Reset function has been activated. The system subsequently switches to passive mode.

## 4.5 Emergency Stop

**! IMPORTANT**  
When the emergency stop is activated, the CPR function can still be activated, by pressing button **(A)**, to lower the wings to horizontal position if required.


When the emergency stop is activated it stops all controls and movements immediately (wings will stay in current position).

1. To activate, push the button **(B)**.
2. To deactivate, rotate the button **(B)** clockwise.
3. Press button **(A)** to reset the system.

## 4.6 Manual operation

**! WARNING!**  
**Risk of squeezing and injury**  
If both wings are simultaneously elevated to angles of more than 12° there is a risk of squeezing the patient.  
– This option must only be used by trained personnel.

**i** The first time MANUAL is activated, after the system has been switched on, it will reset itself before the function can be used.

1. Press and hold button **(H)** for 3 sec. to activate manual operation.
  2. Press buttons on hand control as required.
    - To raise the left wing, press **(I)**
    - To lower the left wing, press **(J)**
    - To raise the right wing, press **(K)**
    - To lower the right wing, press **(L)**
- i** -Left and right buttons can be activated at the same time to move both wings simultaneously.  
-The wings can be moved in a range from 0° to 60°.
3. To deactivate manual operation, press button **(H)** or the ON/OFF & CPR button **(A)**.

### Status indicators

Display	LED
0 = no activity	MANUAL LED active
U = wing moves up	
d = wing moves down	


### 4.6.1 Recommended Procedures

#### Transfer from flat position to the side:

1. Raise both wings simultaneously to increase the comfort for the end user.
2. When both wings are slightly raised, continue raising only the primary wing and leave the secondary wing in position.
3. When the primary wing has reached the desired height, you can (but don't have to) lower the secondary wing to horizontal.

#### Transfer from one side to the other:

1. Simultaneously raise the lower side and lower the higher side
2. Stop lowering the secondary wing when it is still slightly raised and continue raising only the primary wing.
3. When the primary wing has reached the desired height, you can (but don't have to) lower the secondary wing to horizontal.


## 4.7 Automatic operation

The AUTO program is pre-defined and can not be changed.

Program 1 and 2 can be individually defined by the care staff (see chapter *Programming*). Upon delivery the programs are empty and need to be defined before usage.

**i** -If one of the programs is activated, the system will reset to FLAT position before starting the selected program.

-If a new program is selected while another program is active, the system will reset to FLAT position before starting the new program.


1. Press program buttons as required.
  - To start the AUTO Program, press and hold **(M)** for 3 sec.
  - To start Program 1, press and hold **(N)** for 3 sec.
  - To start Program 2, press and hold **(O)** for 3 sec.
2. To stop a running program press the ON/OFF & CPR button **(A)** or press button **(M)** **(N)** or **(O)**, depending on which program is currently active.

**i** The hand control will not work when a program is active.

Status Indicators

Active Program	Display	LED
AUTO	AU	AUTO LED active
Program 1	P1	P1 LED active
Program 2	P2	P2 LED active

4.7.1 Auto program sequence

- Starting from the flat position, both wings immediately begin to raise simultaneously.
- Left wing stops at 12° and the right wing stops at 30°.
- 5 seconds after the right wing has reached 30° the left wing lowers to 0°.
- This position A is maintained for 30 minutes.
- Simultaneously the left wing begins to raise and the right wing begins to lower.
- The right wing stops at 12° and the left wing continues to raise up to 30°.
- 5 seconds after the left wing has reached 30° the right wing lowers to 0°.
- This position B is kept for 30 minutes.

The cycle will continue shifting between position A and B, until the program is stopped.

## 5 Programming

### 5.1 Individually defined programs

Program 1 and 2 can be individually defined by the care staff by setting the angles of the wings and the times for holding the different positions.

#### Defining angles

The angles of the left and right wing in the first side position can be defined individually. The second (opposite) side position will automatically be defined as a mirror image of the first side position.

- If the angle for one wing is set between 13°–30°, this wing will be defined as the “primary wing”.
- The opposite wing will then be defined as the “secondary wing” and limits its selectable angle values to 0°–12°.
- The “primary wing” defines the first side position in the program sequence.

Two options for the wing the end-user is laying on (e.g. left wing if the end-user is laying on the left side):

- If the angle of the secondary wing is set to 0° it will be elevated by 12° during the turn movement and move to a horizontal position 5 sec. after the primary wing has reached its set angle, so the end-user will rest on a flat secondary wing.
- If the angle value for the secondary wing is set between 1°–12° it will stay elevated at the selected angle in a side position and the end-user will rest on a elevated secondary wing.

#### CAUTION!

##### Risk of falling

If the secondary wing stays elevated while the end-user is resting in a side position, the distance to the top of the side rail is decreased.

- Only use the side rails listed for this situation in *9 Technical Data, page 21*.

#### Defining the time:

The times for the two side positions and the flat position can be defined individually between 000 and 180 minutes.

The program sequence always starts in the flat position.

If the time for one side position is set to 000 minutes:

- the program will start to transition from the flat position into the first position after the specified time.
- the program will only cycle between the flat and the other side position (e.g. flat-right-flat-right-...)


If the time for the flat position is set to 000 minutes:

- the program will start immediately after activation with the transition to the first side position.
- the program will only cycle between the two side position (e.g. left-right-left-right-...).

If the time value for the flat position is set between 1–180 minutes:

- the program will start to transition from the flat position into the first position after the specified time.
- the program will always move into a flat position between the two side positions (e.g. flat-right-flat-left-flat-right-...)



## 5.2 How to program the individual automated programs

- Ⓟ Set time / Start programming mode button
- Ⓞ Set angle button
- Ⓡ Adjust button Plus (increase values for time and angle)
- Ⓢ Adjust button Minus (decrease values for time and angle))
- Ⓣ Position indicators (LEDs)
  - Right (top) = right wing or position
  - Flat (middle) = both wings in horizontal position
  - Left (bottom) = left wing or position

First the time during which the wings are in an elevated or flat position for Program 1 or Program 2, is programmed and then the angles.

Step	Function	Display	LED	Note	
1	Press and hold for 3 sec.		SEL	none	Start the programming mode
2	Press		P1	P1	Choose the program to be set
	or		P2	P2	
3	Press		010	P1 or P2 & Right 	The time value appears in the display
4	Press	 	001 020	P1 or P2 & Right 	Adjust the time between 000 and 180 minutes (0, 1, 10-180 in +/- 10 min. intervals)
5	Press		010	P1 or P2 & Flat 	The time value appears in the display
6	Press	 	001 020	P1 or P2 & Flat 	Adjust the time between 000 and 180 minutes (0, 1, 10-180 in +/- 10 min. intervals)

Step	Function	Display	LED	Note
7	Press 	010	P1 or P2 & Left 	The time value appears in the display
8	Press  	001 020	P1 or P2 & Left 	Adjust the time between 000 and 180 minutes (0, 1, 10-180 in +/- 10 min. intervals)
9	Press 	A	None	single *BEEP* = The time has now been set.
10	Press 	015	P1 or P2 & Right 	The angle value appears in the display
11	Press  	014 016	P1 or P2 & Right 	Adjust the angle between 00 and 30 degrees (+/- 1 degree interval)
12	Press 	010	P1 or P2 & Left 	The angle value appears in the display
13	Press  	009 011	P1 or P2 & Left 	Adjust the angle between 00 and 30 degrees (+/- 1 degree interval)
14	Press 	—	None	double *BEEP* = The angle has now been set and the entire program is saved.

-  • The programming mode will automatically be closed down without saving, if no buttons are pressed for 2 min.

### 5.3 Program examples for automated pressure relief

- End-user turning from right side to flat to left side and always laying on a flat wing.
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 0°
  - Time for right, flat and left position > 0 min.
- End-user turning from right to flat to left side and always laying on a slightly elevated wing.
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 1°-12°
  - Time for right, flat and left position > 0 min.

This Program requires “high” side rails.

- End user turning from right to left side, never laying in a flat position and always laying on a flat wing.
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 0°
  - Time flat position = 0 min.
  - Time left and right position > 0 min.
  - see section *Auto program sequence* as an example.

- End user turning from right to left side, never laying in a flat position and always laying on a slightly elevated wing (Cradle function).
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 1°-12°
  - Time flat position = 0 min.
  - Time left and right position > 0 min.

This Program requires “high” side rails.

- End user turning from one side to flat and back to the same side, always laying on a flat wing.
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 0°
  - Time flat position > 0 min.
  - Time left or right side position = 0 min.

- End user turning from one side to flat and back to the same side, always laying on a slightly elevated wing.
  - Angle primary wing = 13°-30°
  - Angle secondary wing = 1°-12°
  - Time flat position > 0 min.
  - Time left or right side position = 0 min.

This Program requires “high” side rails.

- Both wings raised to equal angles without repositioning (Backward hug).
  - Angle primary and secondary wing = 10°-12°
  - Time flat position = 0 min.
  - Time left and right position > 0 min (set values not relevant, see below)

No automatic repositioning will occur when both wings are set to same angle and therefore, the times set for left or right side position do not have any influence.

## 6 Maintenance

### 6.1 General Maintenance Information

- ! **IMPORTANT !**  
Service and maintenance of the Soft Tilt must be carried out together with the bed.
  - For detailed information on maintenance procedures and checklists see *Bed Service Manual* for the bed in use.

### 6.2 Cleaning and disinfection

- ! **IMPORTANT!**  
Wrong fluids or methods can harm or damage your product.
  - Never use corrosive fluids (alkalines, acid, cellulose thinner, acetone etc)
  - Never use a solvent that changes the structure of the plastic or dissolves the attached labels.

#### Electrical components

- ! **IMPORTANT!**  
The IP classification determines the washability of the electrical components.
 

**Electronics classified IPx5 and/or IPx6** may NOT be washed with jet based cleaning equipment or in a washing tunnel.

  - The IP classification is stated on the serial labels of the electrical components.
  - Please also note that the components can have different classification.

### General Cleaning Method

Method: Wipe off with a wet cloth or soft brush.

Max temp: 40 °C

Solvent/chemicals: Mild household detergent or soap and water.


- ! **IMPORTANT**  
For detailed information of cleaning methods see instructions in the section *Non-washable components* in the *Bed User Manual*

### 6.3 Battery

When the Soft Tilt is connected to the power supply, the battery will be fully charged after 12 hours.

We recommend replacement of the battery after 4 years – Call dealer/technician to replace battery.

We recommend a check of the battery every 6 months:

1. Elevate both wings and unplug the power cable.
2. Press ON/OFF & CPR button  and ensure both wings move into horizontal position.

#### Storage

- ! **IMPORTANT**  
Before storage of the Soft Tilt, ensure the system has been turned off completely and the emergency stop button is pressed to prevent the battery from discharging.

---

## 7 After Use

---

### 7.1 Waste disposal

 Waste disposal/recycle must comply with the laws and regulations for waste handling in each country.

Invacare is continuously working towards ensuring that the company's impact on environment, locally and globally, is reduced to a minimum.

We comply with the current environment legislation (e.g. WEEE and RoHS directive).

We only use REACH compliant materials and components.

- All electric parts must be dismantled and be disposed of as electric components.
- Plastic parts must be sent for incineration or recycling.
- Steel parts must be disposed of as waste metals.

**!** **IMPORTANT!**  
• Accumulator back-up  
– Old accumulators must be returned to Invacare or recycled as car batteries

## 8 Troubleshooting

### 8.1 Troubleshooting the electrical system

Error	Display (flashes at 0.5 sec. intervals between two codes)		Acoustic Alarm	LED (flashes red)	Comment	Solution
Stoppage on motor LEFT	E01	E0	10 beeps with 1 sec. ON and 1 sec. OFF  When trying to activate motor again: 1 beep for 1 sec.	Left	only in Auto mode (no display, acoustic alarm and flashing LED in manual mode)	Check the cables to and from the motor for the following: - Cable undamaged. - Cable correctly connected to motor. - Cable connector correctly installed and intact.
Stoppage on motor RIGHT	E02	E0		Right		
Battery mode:  Power failure OR Emergency stop activated	E03	E0	4 beeps with 1 sec. ON and 1 sec. OFF		- 1st acoustic alarm when disconnected from power supply OR the emergency stop is activated. The system switches to battery mode and the wings will stop moving and stay in current position. - 2nd acoustic alarm after 20 min. - 3rd acoustic alarm after 40 min. - After 60 min the system switches to low battery error (E04).	Connect to power supply AND / OR deactivate the emergency stop to charge and check the battery (refer to chapter <i>Maintenance</i> for further information).
Low battery	E04	E0	4 beeps with 1 sec. ON and 1 sec. OFF		- 1st acoustic alarm when the battery is detected to be low while in battery mode. - Further acoustic alarms occur every 20 min until the system turns off to save battery power for one CPR lowering.	
Squeeze detection on motor LEFT	E05	E0	10 beeps with 1 sec. ON and 1 sec. OFF	Auto, P1 or P2 (depending on the interrupted program).	The wing on the side of squeeze detection will elevate to 30° and the opposite wing will move to horizontal position.  - Wing(s) below 12° degree will elevate to 12°. - Wing above 12° will elevate to 30°	Check if any objects are interfering with the wing when moving down and remove.
Squeeze detection on motor RIGHT	E06	E0				
Squeeze detection on motor LEFT and RIGHT	E07	E0				
System can not be turned on					Although the system is connected to power supply it can not be turned on.	Check if the emergency stop is activated and deactivate if necessary.

 Press the ON/OFF & CPR button to reset an ERROR and try again.

Contact your dealer or Invacare representative if the above does not solve your problems.

## 8.2 Programming Errors

Error display for incorrect setup of values for Program 1 or 2

Error	Display	Acoustic Alarm	Comment
Program 1 or 2	E40 for 4 sec.	2 beeps of 5 sec. ON and 1 sec. OFF	<p>Occurs when activating Program 1 or 2 with incorrect time values set.</p> <p>The following combination of time values for the left, right and flat position are defined as incorrect setup:</p> <ul style="list-style-type: none"> <li>- All three time values = 0 min</li> <li>- Only one time value &gt; 0 min / two time values = 0 min</li> </ul>

## 9 Technical Data

### 9.1 Characteristics

Automatic Operation	yes
Max. angle in Manual mode	60°
Max. angle in Automatic mode	primary wing = 30° secondary wing = 12°
Angle values in Automatic mode	0°–30° (+/- 1° interval)
Time values in Automatic mode	0, 1, 10–180 min (+/- 10 min intervals)

### 9.2 Weights

	Max. user weight (provided that the weight of the mattress does not exceed 15 kg)	170 kg
	Max. safe working load weight (including mattress)	185 kg

#### Weights of Soft Tilt components

Complete (main module and 4 arms)	25 kg
Head arm (1 piece)	1.4 kg
Foot arm (1 piece)	2.4 kg
Main module	14.5 kg
Control Box	2.9 kg

### 9.3 Allowed mattress sizes

#### Mattress height and side rail compatibility

**!** **IMPORTANT**  
Depending if the end-user in a side position is laying on a flat wing (mattress horizontal) OR if the wing the end-user is laying on is raised up to 12°, different combinations of mattresses and side rails are required (see appropriate chart below).

#### – End-user laying on a flat wing

Allowed mattress heights, when the end-user is laying on a flat wing (secondary wing = 0°).

Side rail	Mattress height
Britt V / Line / Line Extended	12–14 cm
Britt V / Line / Line Extended with side rail height extender mounted*	15–29 cm
Scala Medium 2	12–16 cm
Scala Decubi 2	17–25 cm

\*Must be dismantled when the side rail is collapsed and the end-user is getting out of the bed.

#### – End-user laying on a elevated wing

Allowed mattress heights, when the end-user is laying on a up to 12° elevated wing (secondary wing = 1–12°).

Side rail	Mattress height
Britt V / Line / Line Extended with side rail height extender mounted*	12–22 cm
Scala Decubi 2	17–18 cm

\*Must be dismantled when the side rail is collapsed and the end-user is getting out of the bed.

#### Mattress width and length depending on bed width

Bed width	Allowed mattress width	Min. length (cm)
85 cm	83–85 cm	200 cm
90 cm	88–90 cm	200 cm

 A castellated mattress is recommended.

### 9.4 Electrical system

Voltage supply: Uin 230 Voltage, AC, 50/60 Hz (AC = Alternating current)
Maximum current input: Iin max.1,5 Ampere
Intermittent (periodic motor operation): Int = Max. 10 %, 2 min ON / 18 min OFF
Insulation class: CLASS II 
Type B Applied Part 
Applied Part complying with the specified requirements for protection against electrical shock according to IEC60601-1. (An applied parts is a part of the medical equipment which is designed to come into physical contact with the patient or parts that are likely to be brought into contact with the patient.)
Sound level: 58.5 dB (A)
Degree of protection: IPx6* or IPx5** The main module, actuators and hand control are protected according to IPx6. The Control box is protected according to IPx5

\* IPX6 classification means that the electrical system is protected against water projected from any direction (not high pressure).

\*\* IPX5 classification means that the electrical system is protected against water low-pressurized water projected from any direction.

## 10 Electromagnetic compatibility (EMC)

### 10.1 General EMC information

Medical Electrical Equipment needs to be installed and used according to the EMC information in this manual.

This product has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2 for Class B equipment.

Portable and mobile RF communications equipment can affect the operation of this product.

Other devices may experience interference from even the low levels of electromagnetic emissions permitted by the above standard. To determine if the emission from this product is causing the interference, run and stop running this product. If the interference with the other device operation stops, then this product is causing the interference. In such rare cases, interference may be reduced or corrected by the following:

- Reposition, relocate, or increase the separation between the devices.

### 10.2 Electromagnetic emission

#### Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group I	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This product is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

### 10.3 Electromagnetic Immunity

#### Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ± 2 kV, ± 4 kV, ± 8 kV, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines; 100 kHz repetition frequency ± 1 kV for input / output lines; 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 0% $U_T$ for 0,5 cycle at 45° steps 0% $U_T$ for 1 cycles 70% $U_T$ for 25 / 30 cycles < 5% $U_T$ for 250 / 300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued operation during power mains interruptions, it is recommended that the product is powered from an un-interruptible power supply or a battery.  $U_T$ is the a. c. mains voltage prior to application of the test level.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 V 150 kHz to 80 Mhz  6 V in ISM & amateur radio bands	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this product is used exceeds the applicable RF compliance level above, this product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating this product.  Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V/m 80 Mhz to 2,7 GHz  385 MHz - 5785 MHz test specifications for immunity to RF wireless communication equipment refer to table 9 of ICE 60601-1-2:2014	Portable and mobile RF communications equipment should be used no closer than 30 cm to any part of this product including cables.

 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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