Andromeda IR-DET

Bruksanvisning / User manual / Brukerveiledning /
Manuel / Gebrauchsanweisung

Art. Nr: 4198XX

Contents

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GB: User manual for Andromeda IR-detectors

Product description

The Andromeda IR-detectors amplify the IR-signals and transfer information to the receiver board. The detectors are available with different casings for in- or outdoor mounting.

419860 Gewa Andromeda IR-DET 3

419865 Gewa Andromeda IR-DET 15

Plastic cover with 3, or 15 m shielded cable.

Discreet mounting that hides the screws.

Surroundings: Indoors and outdoors

Degrees of protection: IP-65

Size: 70 x 42 x 28 mm.

Colour: Polar white.

419870 Gewa Andromeda IR-DET

419872 Gewa Andromeda IR-DET V2

Detector for mounting on a wall.

Surroundings: Indoors.

Degrees of protection: IP-21.

Size: 96 x 96 x 36 mm.

Colour: Polar white.

419875 Gewa Andromeda IR-DET OEM

Detector to be built in.

Size: 49 x 20 x 19 mm.

Cable length: 200 mm.

Material: Metal covered plastic.

Location

If possible, choose a location where the detector is protected against direct sunlight, which reduces the range. Please note that this is an optical system and consequently there must be no obstructions between the transmitter and the detector. Indoors the light however is reflected against walls, floor, and ceiling. Light materials are good reflectors while dark materials are bad ones. Glass will let the IR-light through.

Connection

⚠ CAUTION: Be careful when connecting the detector to the receiver. Making the wrong connection may damage the detector.

NOTE:

- When mounting the Gewa Andromeda IR-DET OEM the detector casing must be isolated from electrically conducting material.
- Detectors included in this manual can only be connected to Gewa Andromeda IR-receivers.

Connect the detector according to the instructions for each individual receiver. The connection instructions refer to the colours of the detector cable. When lengthening the detector cable EKKX 1 x 2 x 0.5 is used. Maximum 100 m detector cable can be connected.

Storage and disposal

Storage

Keep the product indoors at room temperature and away from moist.

Recycling Waste Electrical and Electronic Equipment

The product contains electronic components that must be recycled. When the product is no longer of use, hand it in to a recycling centre for waste electrical and electronic equipment.

Important information

Before you use the product, read and understand the contents of this user manual.

Safety definitions

 **WARNING:** Risk of injury or death if the instructions are not followed.

 **CAUTION:** Risk of damage to the product if the instructions are not followed.

NOTE: Information that is necessary in a given situation.

Serious incidents

If a serious incident related to the product occurs, report this to the manufacturer and the competent authority of your country. A serious incident directly or indirectly led, might have led or might lead to the death of a patient, user or other person; the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

Conformance with requirements for medical devices

The product is marketed as a technical aid for people with disabilities. It complies with all necessary requirements, regulations and directives for medical devices.

Product liability

The product is developed and risk-assessed according to ISO 14971. The manual and the technical specifications are made according to the safety assessments in the risk analysis. Always consider the safety when using the product. Abilia AB does not take responsibility for any consequences following incorrect use or installation of the product. Any form of damage or tampering with the product invalidates the warranty and disclaims the manufacturer's liability. The estimated lifetime of the product applies when used in accordance with its intended use and maintenance instructions.

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Abilia AB, Råsundavägen 6, 169 67 Solna, Sweden
Abilia UK Ltd, Unit 12 Buckingham Business Park, Swavesey,