
DK	Ståsejl og Hoftestøttesejl	2
GB/US	Standing Sling and Hip Support Sling	6
DE	Standstz und Hüftstütze	10
SE	Ståsele och stödsele för höften	14
FR	Harnais de station debout et harnais de maintien des hanches	18
IT	Imbragatura per posizione eretta e imbragatura di supporto dell'anca	22

User manual – vers. 101.00

MD CE

Guldmann™

## DK... STÅSEJL OG HOFTESTØTTESEJL

Vers. 101.00

Varenr.:

Ståsejl  
28336 S/M  
28337 L  
28339 XL  
28340 XXL

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28338 One size

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**GB STANDING SLING AND HIP SUPPORT SLING**

Vers. 101.00

Item nos:

Standing Sling  
 28336 S/M  
 28337 L  
 28339 XL  
 28340 XXL

Hip Support Sling  
 28338 One size

<b>1.00</b>	<b>Purpose and use</b>	<b>6</b>
1.01	Manufacturer	6
1.02	Intended purpose	6
1.03	Area of use	6
1.04	Conditions of use	6
1.05	Important/Precautions	6
1.06	UHF RFID tag	6
1.07	Labels and Marking	6
1.08	Use	7
<b>2.00</b>	<b>Maintenance</b>	<b>7</b>
2.01	Cleaning	7
2.02	The owner's daily maintenance duty	7
2.03	Disposal of slings	7
<b>3.00</b>	<b>Service and lifetime</b>	<b>7</b>
3.01	Safety/service inspections	7
3.02	Lifetime	8
<b>4.00</b>	<b>Technical specifications</b>	<b>8</b>
<b>5.00</b>	<b>EU-Declaration of conformity</b>	<b>8</b>
<b>6.00</b>	<b>Environmental policy statement – V. Guldmann A/S</b>	<b>8</b>
<b>7.00</b>	<b>Warranty and service conditions</b>	<b>9</b>
A.	Warranty	9
B.	Service or Repair	9
	<b>Placing the sling</b>	<b>26</b>

**1.00 Purpose and use****1.01 Manufacturer**

V. Guldmann A/S  
 Graham Bells Vej 21-23A  
 DK-8200 Aarhus N  
 Tel. + 45 8741 3100  
 www.guldmann.com

**1.02 Intended purpose**

The sling is intended for lifting or supporting a person or body parts of a person.

**1.03 Area of use**

The sling is suited for use in hospitals, nursing homes, institutions, rehabilitation centers and in private homes.

The sling is used in combination with Guldmann GLS5 active lifter to assist the user from a sitting to a standing position.

**1.04 Conditions of use**

The sling and the active lifter are designed to lift or transfer a person who is capable of actively participating in the process. The user must also have good communication skills.

The use of the sling is subject to the following:

- The user can balance in a sitting position and has strength in the lower and upper body and can communicate
- The sling is used by trained staff or persons who have been instructed in the use of the sling in question.
- The correct size of sling is used.
- The maximum nominal load of 205 kg must not be exceeded.
- The helper pays attention to the well-being of the user when using the sling.
- The sling is used with the active lifter from Guldmann.

**1.05 Important/Precautions**

- Read the instructions carefully before using the sling.
- The sling's maximum load must never be exceeded.
- The sling may only be used to lift a person.
- Before a sling is used, it must be examined according to point 2.02.
- Never use a sling that is too big for the user.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

**1.06 UHF RFID tag**

This product features a passive UHF RFID tag integrated in the product label. The RFID tag can be used for asset management & tracking purposes. The RFID tag is readable with equipment compliant to EPC global UHF Class I Gen ISO 18000-63.

**1.07 Labels and Marking**



CE marking



Medical Device Class I in accordance with EU MDR Regulation

Read the manual before use

**Example of serial number label****Inspection label**

**Product label****Size label**

 Tumble-drying at low temperature

 Do not iron

**2.02 The owner's daily maintenance duty**

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

**Sling inspection checklist**

Before using a Guldmann sling / accessory check the following:

**Is the sling clean?**

Follow facility specific infection control procedure.

**Is the sling's label present, legible and complete?**

Missing, illegible or incomplete sling label(s) could make identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling impossible.

**Are the lifting straps and stitches intact?**

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

**Is the fabric intact?**

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

**Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?****Conclusion**

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

**1.08 Use**

***If there is any doubt about the selection or use of a lifting sling, please contact your supplier.***

**Important!**


Plan the transfer. Avoid leaving the user in the sling unattended. Do not start to lift until it has been checked that the user cannot get trapped and that the sling does not catch on the bed, wheelchair etc. The user's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes and wires that are attached to the user and/or equipment.

Check that the hand control and hand control cable is free of the lifter, user and other objects before the lifting boom is activated up or down moved.

Guldmann shall not be liable for faults or accidents due to incorrect use of the lifting sling, or for reasons of inadequate attention on the part of the carer or user. If the sling is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

**Caution!**

**Be careful when attaching the sling's lifting straps on the hooks of the lifting boom. Check that the straps have been correctly placed on the hooks. When pressing the up button on the hand control to lift the user, check again that all straps remain correctly placed in the lifting boom's hooks.**

Applying the sling, see page 26.

**2.00 Maintenance****2.01 Cleaning**

Normal washing at the indicated temperature

Do not use bleaching agent

**2.03 Disposal of slings**

Slings are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

**3.00 Service and lifetime****3.01 Safety/service inspections**

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection **must** be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

**Safe Operating Practices with Slings**

Considerations for damaged or defective slings and taking them out of service:

Withdraw the sling from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

**Sling inspection system**

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

**1) Initial**

This level of inspection is done at the time that the sling is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the sling work load limits match those contained in the manufacturer’s catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage.

**2) Frequent**

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.

**3) Periodic**

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling.

Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

**Sling inspection technique**

The sling inspection procedure should be thorough, systematic and consistent; both visual and “hands on” inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the unique sling identification number (important in differentiating similar slings), as well as the condition

of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

**Sample visual examples of synthetic sling damage \*)**


*Chemical/caustic burns*

*Broken stitching*


*Crushed / Frayed webbing*

*Knots*

*Melting / Charring*

*\*) sample visual images not intended to represent all types of potential damage*

**3.02 Lifetime**

The life of the sling is individual and depends on how it is used, washed etc. Before use the sling must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded if necessary.

**4.00 Technical specifications**

Lifting capacity, SWL . . . . . 205 kg (560 lbs)  
 Material . . . . . Polyester

**5.00 EU-Declaration of conformity**

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device class I.

**6.00 Environmental policy statement – V. Guldmann A/S**

Guldmann is continuously working towards ensuring that the company’s impact on the environment, locally and globally, is reduced to a minimum.

It is Guldmann’s goal to:

- Comply with the current environmental legislation (e.g. WEEE and REACH directives)
- Ensure that we, at the widest possible range, use RoHS compliant materials and components
- Ensure that our products do not have an unnecessary negative impact on the environment regarding use, recirculation or disposal
- Ensure that our products contribute to a positive working environment in the places they are utilised

Inspections are made annually by the Department for Nature and Environment from the Municipality of Aarhus using the Danish Environmental Protection Act, section 42 as a reference.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

---

## 7.00 Warranty and service conditions

---

### A. Warranty

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

#### Only for USA

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann's Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

---

### B. Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

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## 4. FR

Approchez le harnais Active Lifter.

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
## 5. GB/US

Place the user's feet onto the foot plate and adjust the knee pad - the top of the knee pad should be level with the top of the shin on a distance of two fingers below the knee joint.

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
## 6. GB/US

**Optional** – if the user has a tendency to have involuntary movements you can use the calf strap to ensure the user's feet do not come off the foot plate. If the calf strap is not in use fasten the strap to the back of the knee pad (avoiding loose straps).

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
## 7. GB/US

Raise the lifting boom to gain tension on the sling – **Ensure that the sling is still secure and attached before continuing with the lift.**

Communication with the user should be maintained throughout the lifting procedure.

Alternatively you have the option to attach the sling onto the lower hooks; this will facilitate an alternative lifting position and is more suitable for the shorter user.

When placing the user back into a seated position push the lifter up to the chair / commode chair / toilet and position the user so the user's calf is parallel to the chair / commode chair / toilet. Begin to lower and ask the user to push his/her pelvis backwards to ensure the correct seated position.

To remove the sling carry out the fitting instructions in reverse order, taking care to remove the sling from the user avoiding skin friction and movement.

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## GB/US

Placing the sling  
– Hip Support Sling with  
the Standing Sling

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Stroppen er nu klar til brug.

## 1. GB/US

**Attachment of the hip support straps**

The support straps have either one or two loops at one end and a buckle at the other. Guide the end with the loop (s) through the hoop of the lifting arm and secure the support strap over the push handle as shown in the illustration.

**Fig. 1 GLS5.1**

On model GLS5.1 hip support straps with 2 loops should be used.

When mounting the hip support straps on the push handle, use the loop located 250 mm from the end as shown in the illustration.

**Fig. 2 GLS5**

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## 4. GB/US

**Optional** – if the user has a tendency to  
have involuntary movements you can use  
the calf strap to ensure the user's feet do  
not come off the foot plate. If the calf strap  
is not in use fasten the strap to the back of  
the knee pad (avoiding loose straps).

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

## 5. GB/US

Attach the hip support sling to the straps snapping the buckle together. Raise the lifting boom to gain tension on both slings – **Ensure that the sling is secure and attached before continuing with the lift.**

Adjust the hip support straps through the buckle and continue. The hip support sling should tension before the standing sling, – if it does not, then lower the user and readjust one or both slings.

As the lift progresses at approximately 60% into the standing procedure the hip support sling will loosen and once in standing the sling will be redundant allowing access to lower garments. Communication should be maintained throughout the lifting procedure.

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## 6. GB/US

When placing the user back into a seated position push the lifter up to the chair/commode chair/toilet and position the user so his/her calf is parallel to the chair/commode chair/toilet.

Begin to lower and at approximately 40% into the lowering procedure re-apply the hip support sling and continue to lower. The sling will ensure correct seated position.

To remove the sling carry out the fitting instructions in reverse order, taking care to peel the sling away from the user avoiding skin friction and movement.

**Alternatively:** The Active Micro Plus sling can be used with the active lifter – for instructions see separate user guide for Active Micro Plus Sling.

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**Guldmann GmbH**
**Guldmann Sverige AB**
**Guldmann Srl**
**Guldmann Inc.**