

ReadEasy *Move 2*

A portable machine that reads to you!

- Easy to use • Clear voices • Transportable
- Accurate • World's fastest reading machine
- Lightweight • Compact • 30 languages



Optional Upgrades: Keypad Feature Pack - Low Vision Touch Pack

VisionAid International

Who is the ReadEasy Move 2 for?

Do you find reading with a video magnifier difficult or impossible? Then ReadEasy Move 2 could be the solution for you.

It translates any printed text into crystal clear speech within seconds. Letters, bank statements, books, recipes, product packaging etc. are all read with ease.

The ReadEasy Move 2 has been specifically designed to be as simple as possible to set up and use. With its intuitive control system comprising just 6 buttons, it is ideal for users of all ages.



Set up & Transport



ReadEasy Move 2 features a unique magnetic quick release camera system. Simply placing the camera into its dock on the base unit automatically locates it into the correct position.

This enables it to be one of the lightest (1.6kg / 3.5lbs) and most compact reading machines available.


Its included padded carry case makes it easy to take with you, wherever you go.

The ReadEasy Move 2 also features an integrated foldable carrying handle making it effortless and safe to move from room-to-room.



Easy to Use

The ReadEasy Move 2 has a specially designed tactile keypad, comprising just 6 buttons making the ReadEasy Move 2 incredibly simple and easy to operate:



Versatility

ReadEasy Move 2's camera allows thick books to be photographed and read accurately, without the need to flatten them, damaging their spines. It also features a special mode for bank statements and tables.



What's more, due to the advanced software used in the Move 2, even round objects such as tins of food can be read aloud.

Use ReadEasy Move 2 for:

- Letters and post
- Bank statements
- Books and magazines
- Newspapers
- Product packaging
- Product instructions
- Recipes and tins of food
- Many more!

Additional Feature Packs



Keypad Feature Pack

For users requiring additional features such as multi-page documents, finer navigation control, saving, loading, importing, exporting of documents & more!

Low Vision Touch Pack

Connecting your Move 2 to a television or computer screen lets you view your documents in large print, (up to 2x newspaper headline size on a 22" screen), as well as listen to them.

The ReadEasy Move 2 can also be used as a magnifier for viewing medicine bottles, photographs & handwritten letters when connected to a television or computer screen.

With a compatible Windows touch screen you can also operate ReadEasy Move 2 with your fingers. Simply pinch to zoom and tap on a word to start reading.



6 Enhanced Text Views

With the Low Vision Touch pack you can view your documents at up to x80 magnification in any one of 6 ways. Text is automatically word wrapped and can be read out at any time.


This makes reading far easier and faster than with a traditional video magnifier and the unique Overlay view ensures even tabular information like bank statements and diagrams are fully accessible.

Image



For photographs and hand written text

Overlay



Ideal for diagrams & bank statements

Column



Most popular view for long text passages

Horizontal




Text on a single line

Vertical



Specifically for tunnel vision (RP)

Word



For larger magnifications

Specifications

Capture Area:	A4, 6 pt text	Power:	110-240 V AC
Camera:	5 Mega Pixel	Connections:	USB 3.0 x 5, HDMI, LAN
Time to read:	2 secs average	Languages:	Auto switching, up to 5
Speakers:	6 Watts RMS	Weight:	1.60 kg (3.5 lbs)
Included accessories:	Headphones, contrast/positioning mat and carry case	Size (WxDxH):	195 x 80 x 215 mm 7.6 x 3.1 x 8.3 inches



VisionAid International Ltd

Bridge Lodge, Spalding Common, Spalding, Lincolnshire, PE11 3AU, England.

Tel: +44 (0)1775 711 977

Web: www.visionaid.com

Email: info@visionaid.com