



## Fitting guide

- ML Vidi is a magnifying galilei system mounted in a frame
  - ML Vidi is easy to adjust to the patient's own prescription and to needs for different viewing distances
  - ML Vidi comes in two different magnifications; 1.6x and 1.8x
  - The ML Vidi system is mounted in a stable frame with a specially made carrier lens
  - ML Vidi is easy to modify with an optional correction lens (for patient's prescription ) behind the system and an optional front lens (for different viewing distances) in front of the system
-

## Fitting procedure for ML Vidi

Follow the steps below in order to get the right solution for your patient.

- Start to put the patient's distance correction in the back of the trial frame
- Put the ML Vidi system in front. Note: The closer the system is to the eye the wider the visual field will be
- With the patient's distance correction placed behind the system, the basic focus distance is 2 m
- Try out the desired magnification. Note: The higher magnification the smaller visual field
- To try out other working distances for the system, use the front lenses available or put additional correction in the back of the system

*Note! It's important that all correction except the front lenses is placed behind the system.*

---

### Example


*Diagnosis: Dry AMD*

*Refraction and acuity: OD:-2.5 -1.5 130 VA: 0.16 Add: 4.0 14p  
OS: -2.0 -1.0 45 VA: 0.2 Add: 4.0 12p*

*This patient found that a magnification of 1.8x was best for TV distance 2 m and for reading a reading cap with power of +4.0 was preferred. Due to photophobia problem ML Filter 511 was preferred.*

## To think about

- To obtain the best possible result, please follow the instructions carefully
- To obtain the best possible visual field, choose a frame where the system can be fitted as close as possible to the eye, we recommend our ML Vidi frame
- To obtain the highest possible magnification, try to reduce the viewing distance as much as possible
- For binocular use, please choose a flat frame for optimal function



- Front lenses with prisms for binocular use is marked for easy identification of prism base

---

### *Ordering of the example on previous page*

*R/L: ML Vidi 1.8*

*Correction:*

*H: -2.5 -1.5 130*

*V: -2.0 -1.0 45*

*ML Filter 511*

*Front lens:*

*R/L: Flip up +4.0 with prisms*

*Mounted in selected frame. Please note the pd and fitting height for correct positioning of the system.*

## Technical data

- Fixed focus galilei-system
- Magnification: 1.6x or 1.8x
- Correction: Sphere -20 to +20 cyl -10
- Default viewing distance: 2 m
- Weight of system without front lenses: 14 gram
- Front lenses available both as fixed model and a flip up model
- The visual field for the two versions are shown below



## Mounting in a frame

- Choose a frame that is stable and can handle the extra weight that comes with the system, e.g. our ML Vidi frame
- We recommend that the order is sent in for mounting at Multilens
- If you choose to do it yourself, see the brochure “Mounting of ML Optical systems”

Multilens AB, Box 220  
435 25 Mölnlycke  
Tel: +46 31 88 75 50  
Email: [info@multilens.com](mailto:info@multilens.com)

