GB Model 600 Installation Manual

DE

Model 600 Installationsanleitung

ES

Modelo 600 Manual de instalación

FR

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IT

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JP

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Stannah

GB Introduction

The information contained within this manual has been designed for use by Stannah approved Technicians / Engineers who have the appropriate product training in the following categories:

1. Product Installation
2. Testing and commissioning

The contents of this manual are in the recommended order of placement, which forms the full installation and commissioning of the Stannah 600.

Pay particular attention to areas marked with this symbol

To reduce the risk of injury, Personal Protective Equipment should be used where appropriate.

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Contents

Installation Tools.....	4
RR Installation Instructions.....	5
Standard Rail installation Instructions.....	7
Setting options.....	13
Chair Installation.....	14
Label Positioning.....	21
Battery charger/remote control: Positioning Guidance Notes.....	22
Installation Overload Test.....	25
Testing and Commissioning.....	31
Customer Handover.....	49
ECU Controls.....	55
Diagnostics.....	56
Adjusting Siena manual swivel.....	74
Rail end - straight finishing dimensions.....	75
Retractable Rail - Adjustment.....	76
Retractable rail – Bus bar cutting Dimensions.....	77
Electrical Diagrams.....	78

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GB

Installation Tools

- 2.5mm 3mm, 4mm, 5mm, 6mm allen key
- Drill driver
- Tape measure
- Extension bar part # 2992266.
- 13mm Ratchet spanner
- 8mm, 10mm Spanner
- Side cutters
- 5 mm, 6 mm Hex key socket
- 13 mm socket
- Torque wrench Range 3Nm—24Nm (2.2lb/ft— 17.6ib/ft)
- PZ3 Bit
- Crimping tool
- Junior Hacksaw
- Spirit level
- Pozi Screwdriver
- Flatblade terminal screwdriver
- Electrical Test meter
- 8mm Twist drill
- 600 rack section

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



Installation Instructions: Retractable Rail

- DE** Montageanleitung: Schiene mit Rückzug.
- ES** Instrucciones de Instalacion: Rail Retráctil.
- FR** Instructions pour l'Installation: Rail Rétractable.
- IT** Istruzioni per l'installazione: Rotaia Retrattile.
- JP** 設置方法: ヒンジレール

- GB**
- Assemble first rail leg to suit hand of rail. Note clamp plate.
 - 1a Fit 1x M4 skt button head screw in position 'A'. Fit rail leg to rail in position arrowed and sit rail on stairs.
 2. Raise retractable arm and fit dynamic section. Torque to **12Nm**
 - 2a Lower Retractable rail into closed position to check latching, using 12v battery
 - 2b Fit transit bracket and tighten.
 3. Place 26mm packer under end of rail. Loosen lower rail leg and allow rail to slide down to packer.
 4. Fit 2nd extended rail leg just below rail joint.
 5. Check rail is upright using spirit level to rear of rail and tighten rail leg to rail.
 6. Fit rail joint pins, rail joint plate and fit and secure top rail section.
 7. Fit uppermost rail leg and Lift top of rail to achieve a clearance of **38mm (1½")** between underside of rail and top nosing.

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
- GB**
8. Check rail dimensions:
 - 8a **170mm (6 3/4")** from front top edge of rail to wall.
 - 8b **38mm (1 1/2")** from underside of rail to stair top nosing.
 9. Fix bottom rail leg to stairs with 1 fixing and remove transit bracket.

Ensure appropriate leg support fixing kits are used for the stair tread type, Eg wood fixings for wooden treads.

10. Fit guide arm to upper Retractable rail arm. Note position!
- 10a Fit slider safety assembly to underside of lower rail section and locate guide arm roller into slot.
- 10b Slide complete assembly down to finish level with motor housing and secure.
11. Power the retractable rail up and down using one 12v battery connected the motor wires to check correct latching. (refer to page 73 if adjustment is required).
12. Fit lower slow/ stop ramp to rear of rail and fit rail end cap.
- 12a Fit lower dynamic cover and refit motor channel cover plate
13. Fit hand-winding bracket to underside of lower rail section close to 2nd rail support.
- 13a Locate hand-winder inside bracket and lower leg support.

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Installation Instructions: Standard Rail

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


- GB 1. Lay lower section of rail on stairs and slide on rail leg with large halfen to bottom of lower rail section

There are two of these legs. The second leg is to be mounted at the top of the upper rail section. Refer to figure 13a.

2. Fit lower slowing / stopping ramp and end cap with spring clip. (Position ramp to rear of rail.) Minimum of 140mm from end of rail.
3. With rail against stair nosings, support rail end 5mm off floor.
4. Check rail is upright using spirit level on rear of rail.
5. Tighten rail leg to rail using 13mm ratchet spanner.

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
GB

6. Fit and tighten 2nd rail leg to rail below joint. Position rail to achieve:

- 7a **150mm (6")** between wall and front of rail.
7b **38mm (1½")** between rail and nosing.
7c **5mm (¼")** between lower end cap and floor.

8. Secure bottom leg with one screw.
9. Fit joint pins and joint plate to lower section of rail.
10. Fit top rail section. Refer to page 75, if cutting on site
11. Fit rail leg with large clamp plate to top of rail for stability.

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
GB

12. Push rail joint together. When satisfied with joint, tighten 3 x joint fixings for security and 2 x cone point grub screws for electrical continuity.
13. Fit further legs as necessary. Ensure:
- 13a **Max 1.5 Metres (59 1/2") between legs. Max 1.0 Metre (39") between rail joints. Minimum 3 Legs.**
- 13b **Fit mini leg at top of rail if supplied, for stability**
14. Screw rail legs to staircase. Ensure rail remains level when securing. If required, fit packers under rail legs to achieve this (part 4303523002 max of 4 per leg).

Ensure appropriate leg support fixings are used for the stair tread type, Eg wood screws for wooden treads. 1.5" / 35mm or 2" / 50mm for Concrete

15. Fit rail leg covers. For Retractable rail, refer to packing list in kit for 1st rail leg cover instructions.

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

GB

16. Remove upper mechanical stop and rack strip stop. Ensure mechanical stop spreader is removed from inside of rail. Remove top section of racking.
 17. Remove carriage cover locking key, carriage covers and upper and lower safety pads.
 - 17a If installing Retractable rail, fit microswitch kit **4301977** to carriage. Refer to fitting instructions within the kit
 18. Using lifting strap, carefully load carriage onto rail.
- Take care not to damage charging brush!**
19. Unclip and lift main ECU. Fit batteries. Switch carriage on via main switch to rear of carriage. Using carriage control buttons, run carriage down rail to a suitable position to fit chair. (Refer to page 53 for control button location.) Switch carriage 'OFF'.

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

GB

20. Refit top section of racking. Ensure pitch is aligned correctly using short section of rack. **Torque to 4Nm (2.95 lb/ft).** Refit mechanical stop and rack strip stop.
21. Fit top slow/stop ramp. Minimum of 140mm from end of rail.
22. Fit battery charger. (Refer to page 22 for location and connection guidance).
23. Slide copper bus bar strip into centre channel at rear of rail. Fit 2 further bus bars if fitting retractable rail. Refer to page 77 for charging and retractable rail bus bar cutting dimensions
24. When bus bar is fully inserted to bottom of rail, slide termination block over bus bar, mark bus bar level with top of termination block and cut.
- 24a Fit self adhesive abrasive pad **4303633** to the middle of the bus bar strip. (Standard rail only).

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- GB**
25. Ensuring terminals are wall facing, connect 24v wire (red) from battery charger.
 - 25a If fitting a retractable rail, feed motor wire to top end of rail under rail legs and secure with spring clips. Cut to length and crimp spade terminals. Connect to terminal block . (Refer to page 79 for connection details)
 26. Slide bus bar back into side extrusion until termination block is positioned against bus bar extrusion. (as photo).
 27. Connect zero volt ring terminal (black wire) from battery charger to underside of rail with shake proof washer.
 28. Cut rack cover to rail length and clip onto rail.
 29. Fit top end cap and secure with spring clip.
 30. If fitting options board, locate and fit 3 spacers and 2 corner fasteners to main carriage ECU tray.
 31. Connect options PCB to main carriage ECU and locate into corner fasteners

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



32

33

- GB
32. Set handing, CON11
Left Hand installation fit link to pins 1&2
Right Hand installation fit link to pins 2&3.
 33. Set carriage options: SW5
A Call chime
B Travel alarm
C Terminus Alarm
D Retractable rail programming switch. (Refer to page 19)
 34. Set chair options PCB (if fitted): SW7
E Seatbelt switch (Isolation).
F Swivel left allowed
G Swivel right allowed
H Footrest request

Viewed from Front

Switches E, F, G are for **Starla** and **Siena** chair options only. For all other chairs, set these switches to the 'OFF' position to ensure correct operation.

SL footrest fold button conversion option only. Set switch "H" to the 'ON' position.

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
- GB**
35. Remove Lower interface from packaging and remove interface covers.
 36. Remove carriage lifting strap and fit bolt in top fixing (arrowed), ensuring spacers are located in correct fixing holes.
 37. Fit lower interface onto top bolt.
 - 37a Fit further two fixings (arrowed), level footrest and tighten interface bolts to **24Nm (17.6 lbf/ft)**.
 38. Fit and secure chair to lower interface. Tighten 4 x fixings to **24Nm (17.6 lb/ft)**

Set seat to footrest height to ensure best fit for user.

39. If setting chair to lowest height position, remove break-off tabs from rear cover.

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
- GB**
40. If fitting manual footrest, remove lower nut from bowden cable adjuster, route bowden cable through lower interface and refit nut.
 41. Feed Bowden cable through cable guide, push cable guide onto footrest pivot shaft and locate cable end in footrest plate.
 - 41a Fit and secure clamp.
 42. Adjust Bowden cable
- Ensure cable guide is located correctly in footrest top plate (arrowed 42)**
43. Fit footrest safety pad. Hook on at top. Grip as fig 43 and press firmly.
 44. Fit footrest carpet

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

GB

Fitting Upholstery

45. **Starla Chair:** Connect seat pressure switch wires and fit seat pad.
- ⚠ 46. Locate front lip and secure to **8Nm (5.9lb/ft.)**.
47. Remove rear casting covers.
48. Fit chair back upholstery and secure to **3Nm (2.2lb/ft.)**.
49. Replace rear casting covers.

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
GB

Fitting Upholstery

1. Fit chair base upholstery (A). Ensure the upholstery is fully pulled over the wings of the seat metal base.
2. Slide upholstery cover (B) over chair back, observing cable route. Fasten the upholstery tabs around the back of the chair.

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
Solus/Sofia Chair Upholstery Fitting

Sofia / Solus

GB

Fitting Upholstery

50. **SL Chair** : Fit chair base upholstery and secure fixings to **10Nm (7.4lb/ft)**.
51. Peel double sided tape off rear of infill pad and fit above seat base. (Sofia only)
52. Remove rear cable access covers and fit chair back upholstery. Secure fixings to **10 Nm (7.4lb/ft)**. Replace rear cable access covers

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


GB

53. Remove 2 x link plugs from carriage ECU.
54. Carefully feed chair cables through carriage aperture.
55. Connect to ECU, referring to the following drawings: **Starla /Siena 4309009 page 81 SL 4309019 page 85**
56. Remove charge brush retaining wire.
57. Programme Retractable Rail if fitted as follows
 - With carriage switched 'ON'
 - Ensure switch 'D' on main PCB is in the 'OFF' position. (refer to page 13 for location).
 - Run carriage up and down over entire rail length. **NOTE!** If powered swivel is fitted, do not allow chair to swivel during this process.
 - ' ' Will be displayed on the seven segment display.
 - Continue to operate until display returns to ' ' or ' '.
58. If required, rehand keyswitch. **Ensure keyswitch is fitted to upside of chair.**
- 58a Connect key switch and refit all carriage and chair covers.

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
GB

59. If Retractable Rail fitted, check clearance between downside safety pad and rail using gauge supplied (59a). If gap is large enough to allow the gauge to slide between rail and safety pad (59b), remove safety pad and fit extension plate to lower safety pad as shown (59c).
60. Adjust upper and lower slowing / stop ramps to desired position.
61. Set measurement at top nosing. Maximum 50mm
62. Apply grease over full length of racking. 20mm at 200mm intervals
63. Fit upper and lower landing remote controls in a safe and convenient position (refer to page 24)
64. Check complete operation.
65. Carry out installation overload test (refer to page 25)
66. Complete Test and Commission sheet, (refer to page 31).
67. Fit carriage labels and if fitted, retractable rail hazard labels. (refer to page 21)
68. Carry out full customer handover, (refer to page 49).

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
Affix Service Branch contact details :

- DE Befestigen Sie die Kontaktinformationen der Zweigstelle
- ES Fije la información de contacto del servicio técnico
- FR Coller les détails de contact du service après-vente
- IT Applicare i dettagli di contatto della filiale di assistenza
- JP サービス支店の連絡先情報を貼付してください


Carriage Label Position: Left Hand

- DE Position der Etiketten an der Fahreinheit: Linksseitig
- ES Posición de las etiquetas del carro: izquierda
- FR Emplacement des étiquettes de chariot : Côté gauche
- IT Posizione dell'etichetta sul carrello: sinistra
- JP 台座ラベルの位置: 左側


Carriage label Position: Right Hand

- DE Position der Etiketten an der Fahreinheit: Rechtsseitig
- ES Posición de las etiquetas del carro: derecha
- FR Emplacement des étiquettes de chariot : Côté droit
- IT Posizione dell'etichetta sul carrello: destra
- JP 台座ラベルの位置: 右側


Affix Manufacturer / Importer label to rear of carriage:

- DE Bringen Sie das Etikett auf der Rückseite des Wagens
- ES Fije la etiqueta del fabricante o de los importadores en la parte trasera del carro
- FR Apposez l'étiquette à l'arrière du chariot.
- IT Applicare l'etichetta sul retro del carrello
- JP キャリッジ後部にラベルを貼ってください


Retractable Rail Hazard Label Position (R/H Shown)

- DE Position der Gefahrenkennzeichnung für die Schiene mit Rückzug auf rechter Seite dargestellt
- ES Posición de la etiqueta de peligro del rail retráctil de la versión de mano derecha
- FR Emplacement de l'étiquette de danger du rail rétractable (côté droit montré)
- IT Posizione dx illustrata dell'etichetta di pericolo della rotaia retrattile
- JP R/H の格納式レール危険標識の位置が示されています


Affix 'Stairlift Supply' label to charger supply socket

- DE Bringen Sie das Etikett „Stairlift Supply“ an der Steckdose des Ladegeräts an
- ES Coloque la etiqueta 'Suministro del salvaescaleras' en la toma de corriente del cargador
- FR Apposez l'étiquette « Stairlift Supply » sur la prise d'alimentation du chargeur
- IT Applicare l'etichetta "Alimentazione montascale" sulla presa di alimentazione del caricatore
- JP 「Stairlift Supply」ラベルを充電器供給ソケットに取り付けます

Public Buildings only: Affix label at each boarding point

- DE Nur öffentliche Gebäude: Etikett an jedem Einstiegspunkt anbringen
- ES Solo edificios públicos: coloque la etiqueta en cada punto de embarque
- FR Bâtiments publics uniquement : apposer une étiquette à chaque point
- IT Solo edifici pubblici: apporre un'etichetta in ogni punto di imbarco
- JP 公共の建物のみ: 各乗り場にラベルを貼る

Battery Charger: Position and Connection Guidance Notes

When positioning the battery charger;

- Connection of the battery charger input cable shall be carried out to the latest edition of IEE wiring regulations or the appropriate national standard,
- Fix the battery charger to a suitable surface using fixing screws provided,
- Mechanically protect battery charger input cable for example, inside trunking,
- **Do not** fix battery charger to stairlift rail,
- **Do not** position battery charger in location which obstructs battery charger ventilation slots,
- **Do not** open battery charger enclosure,
- Any work to remove / change the battery charger input plug (if fitted) must be carried out by an appropriately qualified person and meet the relevant national standard,
- **IMPORTANT NOTE : removing / changing the input plug as supplied may invalidate standard approval of battery charger.**
- Avoid positioning the battery charger in location where susceptible to accidental mechanical damage for example, kicked or struck with vacuum cleaner,
- Avoid positioning the battery charger in location where accidental exposure to fluid is likely for example, ground level where floor is cleaned with fluid,
- Avoid positioning the battery charger where it will cause inconvenience to pedestrian users of the staircase.
- When connecting 24v DC cable, cut to correct length. Crimp connections if required as follows: Red +24v = Spade terminal / Black 0v = Ring terminal

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GB Standard Controls

- Fix holster (2) to wall or newell post at each boarding point
- Fit control (1) into holster
- In order to comply with EN81-40, control must be fixed into holster using screw (3)
- NOTE: For USA/Canada the Controls must be fitted at a height between 890mm—1200mm

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

Overview

It is necessary to complete an important overload safety check upon completion of the installation. A simplified test procedure in accordance with European standard EN81-40 has been created for installed products*.

Details of the test procedure to be done on completion of installation.

1. You should conduct a safe dynamic test by sitting in the chair and riding the stairlift over its full length of travel. Particular attention should be given to rail security and stability.
2. To overload the stairlift a static test should be conducted. You should park the stairlift at the bottom park position, secure your toolbox onto the chair or arm rests, then stand on the front edge of the footrest to create maximum loading, see illustration. Observe the rail security and stability.
3. The same type of static overloading test should be conducted at the upper park position and a mid point on the rail, providing there is a suitable handrail or feature to give you a hand hold while conducting this test at height. Observe the rail security and stability.
4. To carry out a static overload test on a retractable rail, park the stairlift directly over start leg support providing there is a suitable handrail or feature to give you a hand hold. Carry out overloading procedure as above.

What to look for when riding stairlift

The rail legs should remain firmly attached to the floor.

The chair should not tilt excessively.

Expected minimum clearances should be maintained throughout journey.

What to look for during overload test

The rail leg should remain firmly attached to the floor.

No excessive flexing of the rail.

Caution!

- When conducting the overload test ensure the toolbox is secured to the chair to prevent it falling.
- Your toolbox could mark the chair surfaces ensure you adequately protect the chair.
- When conducting the overload test at height ensure you are able to hold something secure to prevent falling.
- When conducting the overload test at height ensure there are no persons directly below you.

*Simplified site test derived from calculations based on rated load and installer weight. (Stannah calculation sheet CS 250-14)

At overload, the product will flex more than usual. The footrest will flex and get closer to the nosings, although it is not expected that it will make contact with the nosings. However, if it does this is acceptable, as in ordinary use the safety pads will ensure that the product stops safely.

A failure is if the overload results in any structural failure of any part of the stairlift or if the fixings pull out, causing a change in the security of any of the rail support fixings. If this occurs, then the stairlift is not suitable for end user use.

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Name and Address of User / Installation site: Tel:

1. Description

a) Manufacturer: Stannah Stairlifts Ltd

b) Stairlift ref number: Chair batch No.: Carriage batch No.: Int/face batch No:

c) Safe Working Load	24 to 52 degrees = 160kg (352lbs)	Siena / Starla Manual and Powered Swivel
	24 to 52 degrees = 135kg (297lbs)	Sofia / Solus Manual and Powered Swivel

d) Rated Speed: 0.1m/s
 e) Motor tested on assembly: YES
 f) Control actuator type: Constant pressure controls
 g) Model type: 600
 h) Type of motor: DC
 i) Motor protection (type): Electronic
 j) Type of drive: Rack and pinion

2. Electrical Tests

a) Confirm that the electrical supply to which charger is to be connected to complies with all applicable local regulations				YES	
b) Battery voltage:		V DC	f) Polarity Check		L-N = L-E= N-E=<30Vac
c) Functional overload protection test:		Yes	g) 0-volt rail connected and continuity tested:		YES
d) Tripping time:		Sec	h) Confirm battery polarity check:		+VE -VE.....
e) Mains fuse rating:		Amps	i) Motor reversal delay:		 Seconds

3. Battery Checks

a) Battery charge voltage (Battery disconnected):	V DC	d) Confirm no leak when battery is fitted:	YES
b) Confirm batteries are Flame retardant and correctly fitted:	YES	e) Confirm off charge indicator is functional:	YES
c) Confirm batteries are fully charge:	YES		

4. Mechanical / Operation Test**SWITCH OPERATION**

a) Confirm chair and wall controls function correctly:		YES
b) Confirm all on / off stop switches operate correctly:		YES
c) Confirm all control limites operate correctly:		YES
d) Confirm final limits operate correctly:		YES
e) Confirm all safety edges operate correctly: As listed below:	Stop	Opposite Direction
e1) Footrest upside safety edge:	YES	YES
e2) Footrest downside safety edge:	YES	YES
e3) Footrest underside safety edge when unfolded:	YES	YES
e4) Footrest underside safety edge when folded:	YES	YES
e5) All carriage upside safety edges:	YES	YES
e6) All carriage downside safety edges:	YES	YES
f) Tests to verify the function of the Overspeed Governor and Safety Gear:		YES
g) Confirm that all interlock switches operate correctly (e.g on rotating seat, downside chair arm isolation switch etc):		YES
h) Confirm retractable rail safety pad switch stops carriage in down direction:		YES
i) Confirm that retractable rail safety arm prevents carriage from clashing with the retractable rail:		YES
j) Confirm that retractable rail clutch slips when manually raised:		YES

COMPONENT FUNCTIONING

k) Ensure that the stopping distance of the carriage with rated load upon application of the brake using a control device is less than 50mm in both directions	YES	
l) Confirm operation of the emergency carriage control buttons	YES	
m) Confirm operation of RR emergency winding function if fitted	YES	
n) Confirm Over Speed Governor has been subjected to factory test for correct tripping speed	YES	
o) Confirm Installation Overload test has been carried out:	YES	
p) If an alarm is fitted, confirm mode of operation and that the system operates correctly:	YES	
q) Confirm clearance of a minimum of 100mm between any part of the stairlift and building structure which is NOT protected by a sensitive edge (safety pad) whilst in travel position.	YES	
r) Confirm clearance of a minimum of 10mm between any part of the stairlift and a continuous smooth surface NOT protected by a sensitive edge whilst in travel position.	YES	
s) Ensure minimum 500mm between adjacent stairlifts during stairlifts during full length of travel:	YES	
t) Confirm that main entrance and exit doors are not obstructed when stairlift not in use:	YES	
u) Confirm user position clearances at any height restriction during full length of travel:	YES	
v) Confirm required hazard warning labels are fitted:	YES	
w) State distance of footrest above floor level: Not to exceed 170mm:	Upper.....mm	Lower.....mm
x) State distance of seat above floor level: Not exceed 650mm:	Upper.....mm	Lower.....mm

5. Confirmation

Confirm operating instructions have been handed to the user and the user has been instructed on use of stairlift: YES / NO

6. Declaration

We certify that on (DATE) This stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.

Signed: Position / Qualifications:

Installation Office / Branch Address:

Introduction

This work method is used for guidance when commissioning a stairlift and completing the Certificate of Test and Examination after the stairlift is installed.

Section 1— Description

- Enter details:
 - Contract or order reference number e.g. 99390
 - The chair serial number.
 - The carriage serial number.
 - The interface serial number.
- Fit correct weight sticker to carriage label from label pack.
- Check that the motor tested label is present on the motor.

Section 2— Electrical Tests

- Confirm that the electrical connection to which the battery charger is to be connected complies with all applicable local regulations, with supporting documentation.
- Measure battery voltage with charger disconnected.
- Carry out functional test to confirm current overload device trips effectively.
- Electronic trip time (maximum 3 seconds).
- Fuse rating (max 2 Amps).
- Polarity check. If a switched, fused spur is fitted, the supply voltage should be measured on L-N, L-E and N-E, using a multimeter set to VAC and the values recorded. If an existing 3 pin socket is to be connected to, a socket tester should be used to validate the socket. In this case, a supply voltage measurement is not required for L-N and L-E, just a tick. An N-E voltage measurement however, is required. On an existing 3 pin socket, any socket tester used must conduct a N-E voltage measurement test and indicate if this is acceptable or faulty. The Martindale EZ165 (others are available), carries out a N-E voltage measurement test and will indicate if the voltage is above 30V AC. If the voltage is above 30V AC the tester will not complete its remaining tests. The 3 pin socket should not be used until an electrician has identified and rectified the fault.
- Carry out a continuity test to ensure the 0 volt wire is mechanically and electrically connected to the rail.
- Check polarity of battery connections.
- Motor reversal delay (min 3 seconds).

Section 3—Battery Checks

- Measure battery charger output voltage on battery wires, with batteries disconnected.
- Batteries are correctly fitted & securely in place & are flame retardant.
- Test battery charge condition using the four volt drop principle.
- Visually check for any leaks or corrosion.
- Check off charge indicator operates.

Section 4 — Mechanical/Operation Test

Switch Operation

- Confirm chair arm and wall controls operate stairlift correctly and in the appropriate direction of travel.
- Confirm all **ON/OFF** switches operate.
- Confirm all stopping limits operate.
- Press ultimate switch on carriage. **9 1-A 5** should show on ECU.
- Ensure the stairlift stops when the safety edges are pressed as listed on the test sheet and can drive away in the opposite direction.

- Test and verify the function of the overspeed governor and safety gear. Refer to page 33 for work method.
- All electrical safety interlock devices function correctly.
- With downside safety pad removed, press retractable rail safety pad switch whilst placing a downward call. Ensure carriage stops.
- With all carriage covers fitted, ensure that the retractable rail safety arm, in the upper, middle & lower position stops the carriage from clashing with the retractable rail. Refer to PDN/ 15/6367.
- Check that retractable rail clutch allows it to be manually raised.

Component Functioning

- Confirm that the stopping distance of the carriage with rated load upon application of the brake, does not exceed 50mm, when using the chair arm controls.
- Check operation of the emergency carriage control buttons.
- Check that the retractable rail can be lifted manually using the emergency winding handle (if applicable)
- Confirm the tested mark is present on the overspeed governor assembly.
- Carry out Installation overload testing. Refer to page 25 for work method.
- Confirm carriage and Retractable rail alarms and warning LED operate (where fitted), including two way voice communication system where applicable.
- Confirm that a minimum *100mm clearance is present between any part of the stairlift and building structure which is NOT protected by a safety edge while in travel position.
- Confirm clearance of a minimum of *10mm between any part of the stairlift and a continuous smooth surface NOT protected by a sensitive edge whilst in travel position.

*Reduced dimensions may be acceptable in certain circumstances where the appropriate documented risk assessments are present.

- Where 2 stairlifts are fitted on the same staircase, ensure that there is a minimum clearance of 500mm between adjacent footrests when unfolded.
- Confirm that main entrance and exit doors are not obstructed when stairlift not in use
- User clearances are measured at = Minimum width of 400mm, measured 200mm in either direction from centre of chair. Minimum depth of 400mm measured from back of chair. Minimum height of 1000mm between seat centre and any height restriction during full length of travel. **Fig 1**
- Confirm that all supplied warning labels are fitted in the correct position, including retractable rail.
- Confirm footrest height at upper and lower floor levels. Record Footrest heights in appropriate box. Maximum footrest height not to exceed 170mm.
- Record seat height. Maximum permissible seat height not to exceed 650mm

Section 5 — Confirmation

The stairlift user must be fully instructed on the use of the stairlift, including the emergency handwinding procedures. All operating instructions should be handed over to the user on completion.

Section 6 — Declaration

This section must be completed in fully using the distributor's address and signed.

On completion of the final test and examination and completion of the test certificate, affix the warning labels to the carriage assembly. This is important to ensure compliance with CE directives.


Fig 1

Tests to verify the function of the Overspeed Governor and Safety Gear

Move the carriage to a safe and convenient position on the rail to carry out the work.

Ensure there are no trapping hazards adjacent to the OSG whilst the carriage is moving and the Safety Pad is removed.

- Using a suitable tool, (suggest screwdriver minimum length 100mm or 4”), gently lift the OSG actuating arm whilst driving the stairlift in an upward direction **via the remote controls at full rated speed** Fig 2, until the OSG switch cuts out the stairlift.
- Check Diagnostics fault code display (**9–1 A-6**) to prove switch disconnection.
- To reset the OSG, hand-wind the stairlift downwards, using a ¼ inch drive socket on the motor hand winding shaft until the OSG cam plate can be seen to spring back into its travel position Fig 1.
- To reset the ECU, switch the rear mounted carriage switch off and then on.
- Drive the stairlift up about 500mm to completely reset the OSG.
- On completion of the work, test the complete operation of the stairlift.


Figure 1

Figure 2

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- 1 Ask user to stand clear of the stairlift making sure that they have a clear view of the lift.
- 2 **Getting to know their stairlift** - Talk the user through the key points of the stairlift .
- 3 **Master power switch / wall switch** - Explain the location and use of these switches. Inform the customer that the stairlift batteries will not charge with the carriage isolator switch in the '0' position.
- 4 **Power key** - Explain the location and use of the power key. Also explain that they have a spare if the original is lost / damaged..
- 5 **Remote calling controls** - Explain the use of the remote calling controls. Also demonstrate how to replace the batteries.
- 6 **Unfolding** - Demonstrate how to unfold the chair and how to avoid the footrest when lowering the seat.
- 7 **Using the stairlift** - Demonstrate how to get in the seat and put the seatbelt on (if immobiliser seatbelt is fitted explain that the lift will not work until this is on). Also explain the need to keep constant pressure in the direction of travel. Demonstrate and explain the arm isolation safety cut out, ensure the user is aware that the downside arm has to be fully lowered in order for the stairlift to operate (if fitted).
- 8 Installers must use a hand held wall control and be in control of the stairlift for the first journey with new user seated correctly in the stairlift chair. Observe and pre-empt any obstructions in the path of the user throughout the journey and take the necessary action to prevent injury to the user. When confident allow and monitor the user operating the stairlift under their own control.
- 9 **Getting out of the stairlift**- Show the user how to operate the swivel levers (or automatic swivel if fitted) and swivel the chair towards the top of the stairs. Ask them to remove the seatbelt and get out of the chair.
- 10 **Folding** - Demonstrate how to fold the stairlift away.
- 11 **Footrest Fold Button** - (if fitted). Demonstrate footrest fold button option.
- 12 Now ask the user to use the stairlift to go down the stairs, covering any of the above that the user seems unsure about.
- 13 Explain the use of the hand-winding controls (including the 2-way automatic swivel if fitted).
- 14 Explain the location of the safety edges and give a demonstration on how they work.
- 15 Check that all of the above has been understood. If they are unsure on any of the above go through it with them again.
- 16 **2-Way Automatic swivel** - (if fitted) Demonstrate how this is used explaining that it will only work when at the bottom of the stairs.
- 17 **2-Way Manual swivel** - (if fitted) Demonstrate how this is used explaining that it will only work when at the bottom of the stairs.
- 18 **Retractable Rail.** Demonstrate how the retractable rail operates, explaining that the area around the retractable rail should be clear before moving the stairlift.
- 19 Explain to the user, that when using the stairlift to ensure any loose clothing is kept clear of any moving parts of the stairlift.
- 20 Explain operation of 2 way communication system.
- 21 Complete handover document.

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
GB

Main ECU controls

- LED9- Status Display** - Operator basic trouble shooting guide
- SW1- Carriage control button up/down**
SW2- Carriage control button up/down
 These buttons allow the carriage to be operated at half speed.
- SW3- Latch switch for normal chair operation and to display user faults (0-9) on fault display LED.**
 De-latch switch to activate carriage control buttons and to enable engineer fault diagnostic mode.
- SW4- Mode Advance button** changes Status display LED (LED9)
- LED's 1- 8 Fault Display.** Illuminate to Identify fault on product.

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
The **Status display** is visible with the carriage cover in place and will show the following basic user fault codes:

- Charge voltage not present.
- 1. Battery voltage too low.
- 2. Ready for operation.
- 3. Safety pad obstructed.
- 4. Chair-arm rocker switch or key switch is OFF.
- 5. Downside chair arm not lowered or seat not loaded.
- 6. Isolator seatbelt not worn. (If Fitted)
- 7. Seat not centred.
- 8. Switch 3 on ECU not latched.
- H-h Switch 3 on ECU not latched. (A stored fault is present, Refer to page 56).
- 9. Fault present. Refer to fault codes.
- F. Software fault. Replace main ECU.


Example Fault:

- ① Stairlift has stopped. Customer reports '9' is displayed in the **Status Display**. *Status Display showing '9' 'Fault present'.*

- ② Remove Carriage Cover to access main ECU and **Fault Display** LED's. Press and de-latch SW 3 on main ECU to enable fault mode.


Status Display continues to show '9' on main ECU.

Status Display	9	A	FAULT
Fault Display LED's	Fault Group 3	1	No power supply detected to option ECU.
		2	Battery voltage in excess of 30V.
		3	Current limit exceeded 24A Instantaneous
		4	Current limit exceeded 16A for 1 second
		5	Powered swivel motor stalled (Starla or Saxon)
			Powered swivel or footrest motor current limit exceeded

Fault Display LED's show fault group '3'. Fault group '3' indicates that it will be one of the faults from 1 to 6.

- ③ For specific fault information, press and hold **Mode Advance** button.

Status display will change from '9' to 'A'.

Status	9	A	FAULT
Fault Display LED's	Fault Group 3	1	No power supply detected to option ECU.
		2	Battery voltage in excess of 30V.
		3	Current limit exceeded 24A Instantaneous
		4	Current limit exceeded 16A for 1 second
			Powered swivel motor stalled (Starla or Saxon)
			Powered swivel or footrest motor current limit exceeded

Fault display changes, to show specific fault. In this case '4'. "Current limit exceeded 16A for 1 second"

Fault display shows specific fault '4'.

If '-', '1' or '2' is shown in the **Status Display** window, removing the carriage cover and delatching SW3 will show **Fault Display LED's** illuminated to indicate which options are fitted to the stairlift.

Fault Display LED's

Status Display '1' or '2'									
LED 1	On	}	On	}	Off	}	LED 5	On	- Power swivel left allowed
LED 2	On		Off		On		LED 6	On	- Power swivel right allowed
LED 3	On - Powered footrest present						LED 7	On = L/H stairlift	Off = R/H stairlift
LED 4	On - Manual 2 way swivel.						LED 8	On	- Option PCB fitted

If '3' is shown in the **Status Display** window, (**Safety pad obstructed**), removing the carriage cover and delatching SW 3 will show **Fault Display LED's** illuminated to indicate which safety pad may be obstructed on the stairlift.

Fault Display LED's

Status Display '3'			
LED 1 On	Bottom footrest safety pad activated		
LED 2	Not Assigned		
LED 3 On	Right footrest/interface safety pad activated		
LED 4 On	Right carriage safety pad activated		
LED 5 On	Bottom interface safety pad activated		
LED 6	Not Assigned		
LED 7 On	Left footrest /interface safety pad activated		
LED 8 On	Left carriage safety pad activated.		

If '9' is displayed in the '**Status Display**' window, refer to the following chart to identify the fault.

Fault Codes

Fault Display LED's

Status Display		FAULT	ACTION
9	A		
Fault Group 1	1	Battery voltage below 23V. Carriage will operate from carriage control buttons only.	Check Charger output / Replace Batteries
	2	Battery voltage below 21V (Carriage will NOT operate from any controls)	Check Charger output / Replace Batteries
	5	Ultimate switch activated	Check for over travel / switch / wiring
	6	OSG switch activated	Check OSG / switch / wiring loom
	7	Ultimate Line switch activated	Check chair arm/ OSG / ultimate switches or looms. Identified ULT 1- 4
Fault Group 2	1	Left and Right slowing switches activated simultaneously	Check switches / wiring loom
	2	Stop Switch activated before slowing switch	Check switch operation / connections
	3	Dual calls from chair arm controls	Check arm control / wiring loom
	4	Dual calls from remote controls	Check IR PCB connections
	5	SL Chair: Seat raised and loaded simultaneously	Check seat switches / loom
	6	'A' Circuit shorted to ground. (Chassis)	Identify and check ' A ' wiring looms
	7	'B' Circuit shorted to ground. (Chassis)	Identify and check ' B ' wiring looms
Fault Group 3	1	No power supply detected to option PCB.	Check fuse F2 (3.15A)
	2	Battery voltage in excess of 30V.	Check charger output / motor
	3	Current limit exceeded 24A Instantaneous	Excess load / Check Motor
	4	Current limit exceeded 16A for 1 second	Excess load / Check Motor
	5	Powered swivel motor stalled	Excess load / Check connections
	6	Powered swivel or footrest motor current limit exceeded	Excess load. Check options motors
Fault Group 4	1	Brake not detected	Check brake and brake connections
	2	Main motor not detected.	Check main motor and connections
	3	Footrest motor not detected	Check footrest motor and connections
	4	Swivel motor 1 not detected	Check swivel motor and connections
	5	Swivel motor 2 not detected	Check swivel motor and connections
	6	Main motor current drops below 1A when travelling upwards	Check motor connections / Check handing plug if travelling downwards

Fault Codes

Fault Display LED's

Status Display	FAULT		ACTION
9	A		
Fault Group 5	1	Starla manual 2 way swivel interlock motor timed out	Check change over mechanism
	2	Starla manual 2 way swivel. Both switches activated Simultaneously.	Check switches / check wiring
	3	Starla / Siena power swivel encoder fault. All outputs 'off'	Check encoder / connections/ slot alignment
	4	Starla / Siena power swivel centred switch fault. Switch still activated when encoder centred	Check switch / encoder connections
	5	Starla power swivel encoder sequence error.	Check / replace encoder
	6	Starla power swivel encoder incorrect sequence in travel	Refer to instruction below to reprogramme encoder / Replace encoder
Fault Group 6	1	Main ECU failure	Replace main ECU
	2	Option PCB failure	Replace option PCB
	3	Intermittent connection between main ECU and option PCB	Check connection between boards
Fault Group 7	1	Con 7 wiring loom short circuit between A&B CCTS	Check wiring loom con 7
	2	Con 9 wiring loom short circuit between A&B CCTS	Check wiring loom con 9
	3	Con 10 wiring loom short circuit between A&B CCTS	Check wiring loom con 10
	4	Con 17 wiring loom short circuit between A&B CCTS	Check wiring loom con 17
	5	Con 18 wiring loom short circuit between A&B CCTS	Check wiring loom con 18
Fault Group 8	1	Chair type Identification different from stored. Reset ECU	Push + Hold SW4 4secs + switch Off/On
	2	Stairlift handing not specified.	Check handing plug fitted correctly
	3	Starla / Saxon swivel wiring incorrect.	Check chair loom con10
	4	New SL swivel detected (SL encoder) not supported.	Replace main ECU to the latest version
F	F	Main ECU software failure	Replace Main ECU
* If Fault Code identifies switch / wiring fault and testing the circuit proves no fault. ECU failure should be assumed.			

Resetting Encoder following a sequence fault

If Resetting the power does not clear an encoder sequence fault, this may be due to the fault occurring during an inwards travel operation. To reset the fault use the following method:

1. Switch power 'OFF'.
2. Set switches 'F' and 'G' on SW7 on options board to the 'OFF' position. (Refer to page 9 for switch location).
3. Switch power 'ON'
4. Switch power 'OFF' and reset switches to desired configuration (Refer to page 9 for switch configuration).
5. Switch power 'ON'.

Accessing Stored Fault Codes

To access this feature:


1. With carriage switched on
2. Remove carriage cover to access the main ECU
3. With status display showing '2', De-latch SW3
4. If there is a stored fault code, the Status Display (LED 9) will cycle between 'H' and 'h'. The Fault Display LED's will cycle simultaneously.

Example: **H4-h3** Indicates '**Footrest motor not detected**' as the last stored fault. (Equivalent to fault 9,4-A,3).

5. To clear the stored fault, press and hold SW4 for 4 seconds. The stored fault will be permanently removed.
6. Re-latch SW3 to return stairlift to normal operation.

NOTE. A stored fault cannot be displayed if the Status Display is indicating '9'.

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
How to reset a Siena manual swivel chair, if it travels in the swivelled position

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
- Loosen the 3 fixings on the switch bracket, as circled above.
- Slide the bracket towards the back of the chair. **NOTE** there is only a small amount of movement.
- Tighten the 3 fixings to 3Nm.
- Test the chair does not run in the swivelled position in all intermediate stop positions.

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Rail end – Straight finishing dimensions

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


Retractable Rail Adjustment

GB Loosen screw 'A' before adjustment of screw 'B'.
After adjustment, re-torque screw 'A' to 12NM.

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






Retractable rail – Bus bar cutting Dimensions

- DE** Schiene mit Rückzug – Schneidmaße der Stromschiene
- ES** Raíl retráctil – Longitud de corte de las barras de distribución
- FR** Rail rétractable – Dimensions de coupe de barre omnibus
- IT** Rotaia retrattile –
- JP** ヒンジレールバス・バー切断寸法



600 RETRACTABLE RAIL CONNECTION DIAGRAM

IMPORTANT NOTES :


- **BATTERY CHARGER MUST NOT BE MOUNTED TO THE RAIL**
- **MAINS INPUT CABLE MUST NOT BE SECURED TO OR STORED INSIDE THE RAIL**

STAIRLIFT SUPPLY WIRING MUST BE CONNECTED IN ACCORDANCE WITH THE APPROPRIATE NATIONAL STANDARD AND THE STANNAH STAIRLIFT INSTALLATION TRAINING PROCESS.


600 DEALER BUS BAR BATTERY CHARGER CONN DIAG

RIGHT HAND

LEFT HAND

600 RETRACTABLE RAIL POWER LOOM ROUTING



FIGURE 1

600 DEALER BUS BAR BATTERY CHARGER & RETRACTABLE RAIL CONN DIAG

FOR DETAILS OF BATTERY CHARGER CONNECTION REFER TO DIAGRAM 4309210.

The diagram illustrates the wiring for the 16-channel digital input module. It shows the following components and their connections:

- KEYSWITCH:** Connected to channels A, B, C, and W. Channel A is connected to PIN 1, B to PIN 2, C to PIN 3, and W to PIN 4.
- BOTTOM FOOTREST:** Connected to channels A, B, C, and W. Channel A is connected to PIN 5, B to PIN 6, C to PIN 7, and W to PIN 8.
- RIGHT FOOTREST:** Connected to channels A, B, C, and W. Channel A is connected to PIN 9, B to PIN 10, C to PIN 11, and W to PIN 12.
- LEFT FOOTREST:** Connected to channels A, B, C, and W. Channel A is connected to PIN 13, B to PIN 14, C to PIN 15, and W to PIN 16.
- RIGHT INTERFACE:** Connected to channels A, B, C, and W. Channel A is connected to PIN 17, B to PIN 18, C to PIN 19, and W to PIN 20.
- LEFT INTERFACE:** Connected to channels A, B, C, and W. Channel A is connected to PIN 21, B to PIN 22, C to PIN 23, and W to PIN 24.
- BOTTOM FOOTREST:** Connected to channels A, B, C, and W. Channel A is connected to PIN 25, B to PIN 26, C to PIN 27, and W to PIN 28.

CHAIR ARM PCB
(MANUAL FOOTREST
VERSION DOES NOT INCLUDE
FOOTREST SWITCH)

The diagram shows a rectangular PCB with a circular footrest switch on the left and four LEDs in the center. Wires connect the switch to a '0' and '1' terminal block. The LEDs are connected to a common ground and individual signal lines.

The wiring diagram shows a 12V battery connected to several components:

- MOTOR 1**: Connected to terminal 10.
- MOTOR 2**: Connected to terminal 9.
- MOTOR 3**: Connected to terminal 8.
- COMMON**: Connected to terminal 12.
- SW A**: Switch A, connected to terminals 1 and 2.
- SW B**: Switch B, connected to terminals 3 and 4.
- GYMR**: GYM relay, connected to terminals 5 and 6.
- BH**: BH relay, connected to terminals 7 and 11.
- LEFT INTERLOCK**: Connected to terminal 1.
- RIGHT INTERLOCK**: Connected to terminal 3.
- INTERLOC**: Connected to terminal 5.
- BE**: BE relay, connected to terminal 10.

[illegible]

Wiring diagram for the LED light fixture. The diagram shows the following components and connections:

- Power Source:** 120VAC (connected to the top terminal of the ULTIMATE LIMIT switch).
- 0VOLT STUD:** Connected to the bottom terminal of the ULTIMATE LIMIT switch.
- Switches:**
 - ULTIMATE LIMIT:** Controls the main power to the light fixture.
 - LEFT SLOW:** Controls the left side of the light fixture.
 - RIGHT SLOW:** Controls the right side of the light fixture.
 - STOP:** Controls the stop function.
- Light Outputs:**
 - RBE (Red Beam):** Connected to the ULTIMATE LIMIT switch.
 - GY (Green Beam):** Connected to the ULTIMATE LIMIT switch.
 - WBR (White Beam):** Connected to the STOP switch.
 - GRN (Green Beam):** Connected to the STOP switch.

LEDs table:

COVER	CARTRIDGE	CARTRIDGE LEFT
1	2	3
2	3	4
3	4	5
4	5	6
5	6	7
6	7	8
7	8	9
8	9	10
9	10	11
10	11	12
11	12	13
12	13	14
13	14	15
14	15	16
15	16	17
16	17	18
17	18	19
18	19	20
19	20	21
20	21	22
21	22	23
22	23	24
23	24	25
24	25	26
25	26	27
26	27	28
27	28	29
28	29	30
29	30	31
30	31	32
31	32	33
32	33	34
33	34	35
34	35	36
35	36	37
36	37	38
37	38	39
38	39	40
39	40	41
40	41	42
41	42	43
42	43	44
43	44	45
44	45	46
45	46	47
46	47	48
47	48	49
48	49	50
49	50	51
50	51	52
51	52	53
52	53	54
53	54	55
54	55	56
55	56	57
56	57	58
57	58	59
58	59	60
59	60	61
60	61	62
61	62	63
62	63	64
63	64	65
64	65	66
65	66	67
66	67	68
67	68	69
68	69	70
69	70	71
70	71	72
71	72	73
72	73	74
73	74	75
74	75	76
75	76	77
76	77	78
77	78	79
78	79	80
79	80	81
80	81	82
81	82	83
82	83	84
83	84	85
84	85	86
85	86	87
86	87	88
87	88	89
88	89	90
89	90	91
90	91	92
91	92	93
92	93	94
93	94	95
94	95	96
95	96	97
96	97	98
97	98	99
98	99	100


**** RIGHT HAND RR CARRIAGE ONLY**

Diagram showing the connection of the 24V output from the battery charger to the 24V output of the centre display.


HEAD OFFICE
Walt Close, East Portway, Andover
Hampshire. SP10 3SD England

Tel : +44(0)1264 364311
Fax : +44(0)1264 339198



INTERFACE & FOOTREST


SIENA WIDE ARM CHAIR

2 WAY POWERED SWIVEL



POWERED FOOTREST




The circuits in this system are in two groups, 'A' and 'B'.
The two groups use a 24VDC pulsing signal.
The 'A' pulses are timed differently to the 'B' pulses.
This adds the electronics detected incorrect wiring and faults.

SL CHAIR




FOOTREST REQUEST SWITCH WIRES
INTO THIS WAGO IN THE CHAIR ARM
SEE DIAGRAM 4309143.






FOR 1 WAY SL. POWERED SWIVEL ONLY CONNECT
1 MOTOR (LEFT OR RIGHT SWIVEL)

INFRARED
TRANSMITTERS

RETRACTABLE



BN MOTOR

SAFETY PAD

LIFELINE

REV	ECN NO.	DESCRIPTION	DATE
0	4247	-	22/03/15
1	4293	RE-SPECIFIED IN STATE 'V' SWITCH - WIRE COILS	26/10/15
2	3569	RE-SAFETY SWITCH DESIGN CHANGED	17/07/13
3	3063	RETRACTABLE RAIL CIRCUIT ADDED	30/08/11
4	3063	FIRST ISSUE (WARNING DESCRIPTION CHANGED)	30/08/11

DRG TITLE : 600 WITH SL CHAIR CONNECTION DIAGRAM				
Stannah Tel : +44(0)1264 364311 Fax : +44(0)1264 339198	HEAD OFFICE		DATE: 07-01-04	
	Watt Close, East Portway, Andover			
	Hampshire, SP10 3SD England			
	CHECKED BY :	DATE:		
	W	IM	KB	DRAWING NUMBER
				REV
			4309019	D


The two groups use a 24VDC pulsing signal. The 'A' pulses are timed differently to the 'B' pulses. This aids the electronics detect incorrect wiring and faults.

RED GRAYE BLACK 24V FUSE BLACK


BRNO COMP TO CHAIRS STUD



CONTROL PCB PULSING CIRCUITS


OPTIONS PCB SIGNAL CIRCUITS

OPTIONS PCB RELAY CIRCUITS

MAIN CONTROL PCB

NO.	SYM.	DESCRIPTION	QTY.
1	RELAY 1	RELAY 1	1
2	RELAY 2	RELAY 2	1
3	RELAY 3	RELAY 3	1
4	RELAY 4	RELAY 4	1
5	RELAY 5	RELAY 5	1
6	RELAY 6	RELAY 6	1
7	RELAY 7	RELAY 7	1
8	RELAY 8	RELAY 8	1
9	RELAY 9	RELAY 9	1
10	RELAY 10	RELAY 10	1
11	RELAY 11	RELAY 11	1
12	RELAY 12	RELAY 12	1
13	RELAY 13	RELAY 13	1
14	RELAY 14	RELAY 14	1
15	RELAY 15	RELAY 15	1
16	RELAY 16	RELAY 16	1
17	RELAY 17	RELAY 17	1
18	RELAY 18	RELAY 18	1
19	RELAY 19	RELAY 19	1
20	RELAY 20	RELAY 20	1
21	RELAY 21	RELAY 21	1
22	RELAY 22	RELAY 22	1
23	RELAY 23	RELAY 23	1
24	RELAY 24	RELAY 24	1
25	RELAY 25	RELAY 25	1
26	RELAY 26	RELAY 26	1
27	RELAY 27	RELAY 27	1
28	RELAY 28	RELAY 28	1
29	RELAY 29	RELAY 29	1
30	RELAY 30	RELAY 30	1
31	RELAY 31	RELAY 31	1
32	RELAY 32	RELAY 32	1
33	RELAY 33	RELAY 33	1
34	RELAY 34	RELAY 34	1
35	RELAY 35	RELAY 35	1
36	RELAY 36	RELAY 36	1
37	RELAY 37	RELAY 37	1
38	RELAY 38	RELAY 38	1
39	RELAY 39	RELAY 39	1
40	RELAY 40	RELAY 40	1
41	RELAY 41	RELAY 41	1
42	RELAY 42	RELAY 42	1
43	RELAY 43	RELAY 43	1
44	RELAY 44	RELAY 44	1
45	RELAY 45	RELAY 45	1
46	RELAY 46	RELAY 46	1
47	RELAY 47	RELAY 47	1
48	RELAY 48	RELAY 48	1
49	RELAY 49	RELAY 49	1
50	RELAY 50	RELAY 50	1
51	RELAY 51	RELAY 51	1
52	RELAY 52	RELAY 52	1
53	RELAY 53	RELAY 53	1
54	RELAY 54	RELAY 54	1
55	RELAY 55	RELAY 55	1
56	RELAY 56	RELAY 56	1
57	RELAY 57	RELAY 57	1
58	RELAY 58	RELAY 58	1
59	RELAY 59	RELAY 59	1
60	RELAY 60	RELAY 60	1
61	RELAY 61	RELAY 61	1
62	RELAY 62	RELAY 62	1
63	RELAY 63	RELAY 63	1
64	RELAY 64	RELAY 64	1
65	RELAY 65	RELAY 65	1
66	RELAY 66	RELAY 66	1
67	RELAY 67	RELAY 67	1
68	RELAY 68	RELAY 68	1
69	RELAY 69	RELAY 69	1
70	RELAY 70	RELAY 70	1
71	RELAY 71	RELAY 71	1
72	RELAY 72	RELAY 72	1
73	RELAY 73	RELAY 73	1
74	RELAY 74	RELAY 74	1
75	RELAY 75	RELAY 75	1
76	RELAY 76	RELAY 76	1
77	RELAY 77	RELAY 77	1
78	RELAY 78	RELAY 78	1
79	RELAY 79	RELAY 79	1
80	RELAY 80	RELAY 80	1
81	RELAY 81	RELAY 81	1
82	RELAY 82	RELAY 82	1
83	RELAY 83	RELAY 83	1
84	RELAY 84	RELAY 84	1
85	RELAY 85	RELAY 85	1
86	RELAY 86	RELAY 86	1
87	RELAY 87	RELAY 87	1
88	RELAY 88	RELAY 88	1
89	RELAY 89	RELAY 89	1
90	RELAY 90	RELAY 90	1
91	RELAY 91	RELAY 91	1
92	RELAY 92	RELAY 92	1
93	RELAY 93	RELAY 93	1
94	RELAY 94	RELAY 94	1
95	RELAY 95	RELAY 95	1
96	RELAY 96	RELAY 96	1
97	RELAY 97	RELAY 97	1
98	RELAY 98	RELAY 98	1
99	RELAY 99	RELAY 99	1
100	RELAY 100	RELAY 100	1

KEY
 (M) = MAIN CONTROL PCB
 (O) = OPTIONS PCB
 (C) = CONNECTOR
 (S) = SWITCH
 (R) = RELAY

DRG TITLE : 600 WITH SL CHAIR KEY WIRING DIAGRAM

HEAD OFFICE
 Wait Close, East Portway, Andover
 Hampshire, SP10 3SD England
 Tel : +44(0)1264 364311
 Fax : +44(0)1264 339198

DRAWN BY :	DATE :
M	03/08/11
CHECKED BY :	DATE :
IW	
DRAWING NUMBER	REV
4309022	C

600617000011

*

'Original instructions in English'