Gowing²

Quick Start Guide / Snelstartgids

ID manual: 606103
Revision / Revisie: 1.0.0
Date / Datum: 01-02-2020

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This is the Quick start guide for Gowing², a dynamic arm support system. This product is developed, manufactured and possibly distributed by Focal Meditech B.V. This quick start guide contains the basic information regarding Gowing². The aim of this information is to ensure successful, safe and effective use of the device. This quick start guide will give you a quick look on the features and important safety warnings of Gowing². Please take note of the safety warnings before using Gowing².

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Overview / Overzicht


Figure 1: Overview Gowling² with terms of specific parts

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Warning: Do not modify any part of this equipment without authorization of the manufacturer. Failure of this can lead to malfunctioning and will void the warranty.

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Warning: In case of a faulty device contact Focal. Do not try to fix it yourself. Failure of this will void the warranty.

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Warning: In case of doubt about the safety of the device contact Focal Meditech.

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Warning: In case of a serious incident when using the device, contact Focal Meditech and the national authority of your country.

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Switch on Gowing² / Aanzetten Gowing²


Figure 2: Picture of on/off switch with indication LED

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Before switching on Gowing², the Gowing² has to be attached to the Power Module (See “Attach / Detach Gowing²”). Gowing² can be switched on by turning the on/off switch to (I). The indication LED will light up green. If the LED lights up red Gowing² might not be properly attached to the Power Module.

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Detach Gowing² from Power Module (right handed version):

STEP 1: Set all settings to default

STEP 2: Remove the user's arm

STEP 3: Switch off Gowing²

STEP 4: Disconnect the human interface and store it in the case

STEP 5: Place your right hand on position 1 and your left hand on position 2 (see [Figure 3](#))

STEP 6: Rotate the red handle upwards in the direction of the arrow (see [Figure 4](#))

STEP 7: Tilt the Gowing² approximately 15° in the direction of the arrow (see [Figure 5](#))

STEP 8: Move the Gowing² up in a slanted line until it is free from the Power Module (see [Figure 6](#))

STEP 9: Store Gowing² always in its case to prevent the Gowing² from breaking or being damaged.

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Attach Gowing² to the Power Module (right handed version):

- STEP 1:** Grab Gowing² with your right and on position 1 and your left hand on position 2
- STEP 2:** Hold Gowing² above the Power Module in a tilted position of approximately 15°
- STEP 3:** Lower Gowing² in a slanted line until it catches the hook of the power module
- STEP 4:** Tilt Gowing² straight until the body of Gowing² is vertical
- STEP 5:** Push down the red handle
- STEP 6:** Attach the human interface to Gowing²
- STEP 7:** Position the user's arm in the arm fitting
- STEP 8:** Switch on Gowing²

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The attach/detach procedure of the left handed version of Gowing² is similar. The only difference is that your left hand is on position 1 and your right hand on position 2.

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The human interface containing the load arm, arm fitting, elbow fitting and wrist support can be removed easily. To detach the human interface, press both red eject knobs at the top of Gowing², and move the human interface away from the body in the upwards direction.

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Attaching the human interface can be done by pushing the lever into the body. The buttons do not have to be pushed. When the human interface is in position, it cannot be removed without pressing the buttons.

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The left column contains the blocking function and the adjustment of the balancing force. Gowing² has blocking mechanisms for its up/down movement, forward/backward movement and rotation (in the horizontal plane). Two buttons are assigned to adjust the balance force in the vertical plane. Pushing the plus (+) button will increase the balance force. Pushing the minus (-) button will decrease the balance force.

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
Figure 9: Adjust tilt angle of Gowling²

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The buttons in middle column are used for the tilt function. The button marked with  will adjust the tilt module to the position where the rotating axis is set perpendicular to the gravity plane. After pressing this button the LEDs around this button will start blinking green and the tilt module will level to 0°. For most users the 0° position will be the best setting. The buttons  for forward and  for backwards can be used to adjust the tilt of Gowling².

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The buttons in the right column are used for the lift function of Gowling². When pushing the button , the human interface will go to the lowest position. The user is now able to move the arm freely again. The buttons  and  can be used to adjust the lift position.

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Eneaswitch² Single switch control / Eneaswitch² een-knops bediening

The single switch mode of Gowing² allows the user to control the Gowing² with only one switch. To control Gowing² in single switch mode there must one external switch be configured as “Single switch”. This can be only be done by authorised people.

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The single switch control starts scanning after pressing the external switch. The function that is indicated by the blue LED will be activated when the external switch is being pressed. For example, In **Figure 11** the “lock” function will be activated when the external switch is being pressed. If the user does not give an input within 1 second the scanning resumes.

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
Figure 12: Lift function selected

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If a function is indicated by two blue LEDs, the corresponding function will also be activated when the external switch is being pressed. Subsequently the activated function will toggle between the two indicated blue LEDs. See **Figure 12** where the lift function will be activated when the external switch is being pressed. If the user does not give an input for 1 second the scanning resumes.

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The complete user manual is available on our website and can be downloaded (<https://www.focalmeditech.nl/nl/documentatie>). In case you want to receive a printed user manual in your own language, please contact us at +31 (0)13 533 31 03 and we will send you without any costs a hard copy.

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