

## Order Form 2017

**Interactive  
with Adobe Reader**


Please refer to the help in the traders area

### Basic model

| X                     | Item no. | Model                                 | Quantity |
|-----------------------|----------|---------------------------------------|----------|
| <input type="radio"/> | 5582-... | <b>Basic model</b><br>Padded strap    |          |
| <input type="radio"/> | 5589-... | <b>Closed top</b>                     |          |
| <input type="radio"/> | 5584-... | <b>With ear protection</b>            |          |
| <input type="radio"/> | 5591-... | <b>Closed top with ear protection</b> |          |

When ordering please forward "Starlight dimension sheet" fully completed.

### Please select your colour

|                                                                                    |                                                                                    |                                                                                    |                                                                                    |                                                                                     |
|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <input type="radio"/> brown                                                        | <input type="radio"/> red                                                          | <input type="radio"/> blue                                                         | <input type="radio"/> green                                                        | <input type="radio"/> black                                                         |
|  |  |  |  |  |



Starlight®  
Flex basic model

### Size selection (please tick as appropriate)

|                                |                                |                                |                                 |
|--------------------------------|--------------------------------|--------------------------------|---------------------------------|
| <input type="radio"/> <b>S</b> | <input type="radio"/> <b>M</b> | <input type="radio"/> <b>L</b> | <input type="radio"/> <b>XL</b> |
| 48 - 52 cm                     | 52 - 56 cm                     | 56 - 60 cm                     | 60 - 63 cm                      |

### Accessories/Options

| X                     | Item no. | Accessories                                                                          | Quantity |
|-----------------------|----------|--------------------------------------------------------------------------------------|----------|
| <input type="radio"/> | 5554-01  | <b>Velcro adjustment instead of Ratchet adjustment, incl. alternative upper side</b> |          |
| <input type="radio"/> | 5555-01  | <b>Chin protection</b>                                                               |          |
| <input type="radio"/> | 5556-01  | <b>Velcro forehead protection</b>                                                    |          |
| <input type="radio"/> | 5557-01  | <b>Velcro sunshield</b>                                                              |          |
| <input type="radio"/> |          | <b>Net surcharge for special colour</b>                                              |          |
| <input type="radio"/> |          |                                                                                      |          |

Additional accessories on next page →

Individual agreements are always superior to standard conditions.

We would like to point out, that we reserve the right to make alterations if required due to technical reasons. Possible alterations made apply also to colours and fabrics. The products displayed are either the basic models or a variation thereof (e.g. with accessories). Please, take this into consideration when ordering. We do not accept liability for printing errors.

The General Terms and Conditions of ATO FORM GmbH apply.

|                                   |       |
|-----------------------------------|-------|
| Address/Stamp or Customer number: |       |
| Name/Company                      | ..... |
| Address                           | ..... |
|                                   | ..... |
|                                   | ..... |
| Customer number                   | ..... |
| Date/Signature                    | ..... |




## Order Form 2017

### Interactive with Adobe Reader

Please refer to the help in the traders area

#### Fasteners

| <input checked="" type="checkbox"/> Item no. | Type                                         | Quantity     |
|----------------------------------------------|----------------------------------------------|--------------|
| <input type="checkbox"/> Standard            | <b>Thorn buckle</b>                          | <b>incl.</b> |
| <input type="checkbox"/> 5182-02             | <b>Fixlock fastener</b>                      |              |
| <input type="checkbox"/> 5182-02-SV          | <b>Fixlock fastener<br/>with safety lock</b> |              |
| <input type="checkbox"/> 5182-02             | <b>Velcro fastener</b>                       |              |
| <input type="checkbox"/> 5182-03             | <b>Magnetic fastener<br/>(including key)</b> |              |
| <input type="checkbox"/> 5182-03-ES          | <b>Spare key</b>                             |              |
| <input type="checkbox"/> 5182-06             | <b>Strengthened thorn buckle</b>             |              |



**The best head protection is the one which is gladly worn. Combined colours, special fabrics and sizes are available on request, for a small additional surcharge.**

Individual agreements are always superior to standard conditions.

We would like to point out, that we reserve the right to make alterations if required due to technical reasons. Possible alterations made apply also to colours and fabrics. The products displayed are either the basic models or a variation thereof (e.g. with accessories). Please, take this into consideration when ordering. We do not accept liability for printing errors.

The General Terms and Conditions of ATO FORM GmbH apply.

Address/Stamp or Customer number:

Name/Company .....

Address .....

Customer number .....

Date/Signature .....

Commission ..... Disease profile .....

**Frequency of falls** ☐ infrequent ☐ often ☐ very frequent


**Fall direction** ☐ forward ☐ backward ☐ left ☐ right

☐ Is there risk of head injury or self-harming behaviour?

Age:


Level of activity:

The head dimensions are determined by using a flexible measuring tape according to the illustrations below.



## General measurements

Fig. 1



\*Occipital protuberance

**A** \_\_\_\_\_ cm  
(Head circumference at eyebrow height, above the occipital protuberance)

**B** \_\_\_\_\_ cm  
(Midway between the eyebrows to upper edge of the occipital protuberance)

**C** \_\_\_\_\_ cm  
(Over the highest point of the head, beginning and ending at the line of head-circumference measurement A)

**D** \_\_\_\_\_ cm  
(Under the chin, beginning and ending at the line of head-circumference measurement A)

Address/Stamp or Customer number:

Name/Company


Address

Customer number


Date/Signature

## Position of the ears on the head


Fig. 2



☐ normal



☐ low



☐ high



In the case of anatomical abnormalities, describe the shape of the head as precisely as possible. Photos are very helpful here.

### Attention


Ideally, the tips of the ears are located on the line of measurement A.

If the positioning of the ears differs from this, make certain that measurement C is taken across the top of the head, beginning and ending at the line of measurement A.

In the case of asymmetrical ear positions, please give a clear and detailed description in the box labelled "comments" below.

## Skull perspectives

Fig. 3



Comments:

If necessary, please print and alter the forehead and/or skull shape in the adjacent illustrations or, if applicable, note any other deviations.

Further particulars, if Varia, Flex, Protect, Protect Plus, Go Sport or Go Sport Plus models are required.

Line of measurement A to lower edge of ears, position of the ears





Fig. 4

**E** \_\_\_\_\_ cm  
(Eyebrow height (line of measurement A to lower edge of earlobes))

**F** \_\_\_\_\_ cm  
(Distance between outer corner of the eye and measurement line E)

Ear measurement

Fig. 5



**G** \_\_\_\_\_ cm  
(Height of the ears)


**H** \_\_\_\_\_ cm  
(Width of the ears)

To be completed if a cheek protector is required!

Size and position of the cheek protector

(For Protect, Protect Plus, Go Sport, Go Sport Plus models)

Fig. 6



**L** \_\_\_\_\_ cm  
(Height of the free area from forehead to cheek protection)


**M** \_\_\_\_\_ cm  
(Height of the cheek protection)

**N** \_\_\_\_\_ cm  
(Width of the cheek protection)

Additional measurements for supply of a full protection helmet

Size of the facial area to be left open

Fig. 7



**O** \_\_\_\_\_ cm  
(Width of the facial area at eye height)

**P** \_\_\_\_\_ cm  
(Height of the facial area from forehead to below the mouth)

Address/Stamp or Customer number:

Name/Company

Address

Customer number

Date/Signature