

General Purpose Series

Our general purpose slings are available in models with head support or without. Both can be used in most lifting situations, fit most patients and are an optimal choice when it is unknown who will be using the sling. Our general purpose slings come in a variety of fabrics, seating positions and sizes to accommodate most needs. The Polyester Net and Trevira options allow the sling to be left behind the patient for a short period of time after the transfer is complete which makes these slings the ideal choice when transferring into a chair or other seating.

Without extra head support



JS JM S M L XL

Basic Sling Model #25045, 25050, 25070

- Available in Polyester Net (25045), Trevira (25050) or Polyester (25070)
- Adapts to the patient's body and offers good support for the entire back
- Trevira and Polyester Net can remain behind patient
- Upright seating position



Model #25045
Polyester Net



Model #25070
Polyester

With extra head support



S M L XL XXL XXXL

Full Body Sling Model #25005, 25010

- Extra head support and divided legs
- Suitable for patients with poor trunk stability
- No handles, reinforcements or padding
- Reclined sitting position
- Model #25005 available only in sizes S to XL




Model #25005
Polyester Net



JS JM S M L XL

Sit Sling Model #25030

- Reinforced back and legs
- Suitable for patients with poor trunk stability
- Optional waist belt
- Divided legs and shoulder height
- Polyester fabric makes sling application fast and simple



JS JM S M L XL

Silhouette Sling Model #25015, 25020

- Extra head support and divided legs
- Narrow overall width make it an ideal option for users who are smaller in stature
- No handles, reinforcements or padding
- Can be left behind patient after transfer is done



Model #25020
Trevira

Sling Solutions

Human Care offers a broad range of anatomically designed slings developed with the safety and comfort of the patient in mind. Our portfolio includes more than 29 different models of slings offered in a variety of sizes and fabrics to accommodate a vast number of applications and needs.

Suitable for:



Choose the right sling

What is the sling being used for?

To/From a bed	Bathing/Toileting	To/From a chair	Gait training/mobilization	To/From the floor	Dressing/Undressing

What size do you need?

Junior - Approx. Patient Weight	Adult - Approx. Patient Weight	Bariatric - Approx. Patient Weight
JS 33-44lb/ 15-20kg	S 66-132lb/ 30-60kg	XL 242-309lb/ 110-140kg
JM 44-77lb/ 20-35kg	M 121-198lb/ 55-90kg	XXL 275-500lb/ 125-227kg
JL 77-110lb/ 35-50kg	L 187-265lb/ 85-120kg	XXXL 350-1000lb/ 159-454kg

Choose the right fabric/finish

Polyester <i>Durable and easy to maintain</i>	Polyester Net <i>Flexible and airy; dries quickly</i>	Trevira <i>Soft and comfortable; flame retardant</i>
Polyester <i>Ventilated for maximum breathability</i>	Cotton <i>Ideal for single patient use</i>	Clip or Loop Sling <i>Two attachment styles</i>

Assess the user

Measurements	User Weight	Assess capabilities
	250kg/551 lb	
Measure the back width and height, throcater to throcater (while seated).	Ensure that the patient does not exceed the maximum weight capacity of these slings. ²	Does the user need assistance? Do they have muscle strength and tone? Are they anxious or have fragile skin? ¹

¹Sampling of questions. See our Sling Evaluation form to assist in assessment of the user's capabilities.
² 25010-XXL, 25010-XXXL maximum user weight 661lb/300kg

Also read our other sling solution brochures: Hygiene Series, Specific Needs Series, Repositioning Series and Bariatric Series.

“ We provide sling solutions to reduce injuries and make safe patient handling possible. ”

Technical Specifications

Basic Sling 25045/50/70*		Part no	A	B	C
Junior Small		25045/*-JS	22cm/8.7"	62cm/24.4"	50cm/19.7"
Junior Medium		25045/*-JM	28cm/11"	65cm/25.6"	54cm/21.3"
Small		25045/*-S	35cm/13.8"	77cm/30.3"	58cm/22.8"
Medium		25045/*-M	40cm/15.7"	88cm/34.6"	63cm/24.8"
Large		25045/*-L	44cm/17.3"	98cm/38.6"	65cm/25.6"
XLarge		25045/*-XL	48cm/18.9"	102cm/40.2"	70cm/27.6"

Sit Sling

Junior Small		25030-JS	22cm/8.7"	56cm/22"	53cm/20.9"
Junior Medium		25030-JM	28cm/11"	63cm/24.8"	58cm/22.8"
Small		25030-S	35cm/13.8"	76cm/29.9"	62cm/24.4"
Medium		25030-M	40cm/15.7"	86cm/33.9"	62cm/24.4"
Large		25030-L	44cm/17.3"	98cm/38.6"	63cm/24.8"
XLarge		25030-XL	48cm/18.9"	102cm/40.2"	64cm/25.2"

Full Body Sling 25005/10

Small		25005-S	26cm/10.2"	76cm/29.9"	86cm/33.9"
Medium		25005-M	26cm/10.2"	85cm/33.5"	95cm/37.4"
Large		25005-L	34cm/13.4"	96cm/37.8"	95cm/37.4"
Extra Large		25005-XL	36cm/14.2"	105cm/41.3"	106cm/41.7"
Small		25010-S	31cm/12.2"	85cm/33.5"	90cm/35.4"
Medium		25010-M	40cm/15.8"	92cm/36.2"	91cm/35.8"
Large	25010-L	40cm/15.8"	98cm/38.6"	100cm/39.4"	
XLarge	25010-XL	44cm/17.3"	105cm/41.3"	102cm/40.2"	
XXLarge	25010-XXL	44cm/17.3"	111cm/43.7"	106cm/41.7"	
XXXLarge	25010-XXXL	44cm/17.3"	118cm/46.5"	108cm/42.5"	

Silhouette Sling 25015-20

Junior Small		25015/20-JS	18cm/7.1"	62cm/24.4"	64cm/25.2"
Junior Medium		25015/20-JM	22cm/8.7"	67cm/26.4"	71cm/28"
Small		25015/20-S	26cm/10.2"	77cm/30.3"	79cm/31.1"
Medium		25015/20-M	29cm/11.4"	87cm/34.2"	90cm/35.4"
Large		25015/20-L	35cm/13.8"	93cm/36.6"	96cm/37.8"
XLarge		25015/20-XL	40cm/15.8"	97cm/38.2"	106cm/41.8"

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone: +46 8 665 35 00
Fax: +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Phone: +61 1300 132 736
Fax: +61 03 9645 4954
info.au@humancaregroup.com

HUMAN CARE NETHERLANDS

Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE CANADA

Phone: 613.723.6734
Fax: 613.723.1058
info.ca@humancaregroup.com

HUMAN CARE UNITED STATES

Phone: 512.476.7199
Fax: 512.476.7190
info.us@humancaregroup.com



HUMAN CARE
Your Life. Your Way.

Sling Solutions
General Purpose Series